

GW01 GLASS BEAKERS & FLASKS

Volumetric Flask With I/C Joint Stopper

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

9
PRODUCTS

9
PACK
OPTIONS

₹650
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Volumetric Flask With I/C Joint Stopper			
GLT-GW01-001	Volumetric Flask With I/C Joint Stopper — Capacity: 5ml	Pack of 2	650
GLT-GW01-002	Volumetric Flask With I/C Joint Stopper — Capacity: 10ml	Pack of 2	800
GLT-GW01-003	Volumetric Flask With I/C Joint Stopper — Capacity: 25ml	Pack of 2	900
GLT-GW01-004	Volumetric Flask With I/C Joint Stopper — Capacity: 50ml	Pack of 2	1,050
GLT-GW01-005	Volumetric Flask With I/C Joint Stopper — Capacity: 100ml	Pack of 2	1,100
GLT-GW01-006	Volumetric Flask With I/C Joint Stopper — Capacity: 200ml	Pack of 2	1,200
GLT-GW01-007	Volumetric Flask With I/C Joint Stopper — Capacity: 250ml	Pack of 2	1,300
GLT-GW01-008	Volumetric Flask With I/C Joint Stopper — Capacity: 500ml	Pack of 2	1,400
GLT-GW01-009	Volumetric Flask With I/C Joint Stopper — Capacity: 1000ml	Pack of 2	1,550

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-01

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Volumetric Flask With I/C Joint Stopper, Evaporating Flask & more

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

8
PRODUCTS

8
PACK
OPTIONS

₹275
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Volumetric Flask With I/C Joint Stopper			
GLT-GW01-010	Volumetric Flask With I/C Joint Stopper — Capacity: 2000ml	Pack of 2	3,200
Evaporating Flask (rotary evaporator type)			
GLT-GW01-011	Evaporating Flask (rotary evaporator type) — Joint: 100ml	Pack of 6	3,900
GLT-GW01-012	Evaporating Flask (rotary evaporator type) — Joint: 300ml	Pack of 3	2,550
GLT-GW01-013	Evaporating Flask (rotary evaporator type) — Joint: 600ml	Pack of 2	2,550
GLT-GW01-014	Evaporating Flask (rotary evaporator type) — Joint: 800ml	Pack of 2	3,900
R.B Flask Single Neck With Joint			
GLT-GW01-015	R.B Flask Single Neck With Joint — Capacity: 25ml · Joint: B-14	1 pc	275
GLT-GW01-016	R.B Flask Single Neck With Joint — Capacity: 50ml · Joint: B-19	1 pc	325
GLT-GW01-017	R.B Flask Single Neck With Joint — Capacity: 100ml · Joint: B-24	1 pc	350

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-02

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

R.B Flask Single Neck With Joint & R.B Flask Three Neck With Joint

Illustrative image

Low - and tall - form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

9
PRODUCTS

9
PACK
OPTIONS

₹300
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
R.B Flask Single Neck With Joint			
GLT-GW01-018	R.B Flask Single Neck With Joint — Capacity: 250ml · Joint: B-24	1 pc	400
GLT-GW01-019	R.B Flask Single Neck With Joint — Capacity: 500ml · Joint: B-24	1 pc	450
GLT-GW01-020	R.B Flask Single Neck With Joint — Capacity: 1000ml · Joint: B-24	1 pc	525
GLT-GW01-021	R.B Flask Single Neck With Joint — Capacity: 2000ml · Joint: B-34	1 pc	1,325
GLT-GW01-022	R.B Flask Single Neck With Joint — Capacity: 3000ml · Joint: B-34	1 pc	3,350
GLT-GW01-023	R.B Flask Single Neck With Joint — Capacity: 5000ml · Joint: B-34	1 pc	4,525
GLT-GW01-024	R.B Flask Single Neck With Joint — Capacity: 10000ml · Joint: B-45	1 pc	9,000
GLT-GW01-025	R.B Flask Single Neck With Joint — Capacity: 20000ml · Joint: B-55	1 pc	11,150
R.B Flask Three Neck With Joint			
GLT-GW01-026	R.B Flask Three Neck With Joint — Capacity: 25ml · Joint: B-14, B-14, B-14 (centre, side necks)	1 pc	300

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-03

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

R.B Flask Three Neck With Joint

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

9
PRODUCTS

9
PACK
OPTIONS

₹400
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
R.B Flask Three Neck With Joint			
GLT-GW01-027	R.B Flask Three Neck With Joint — Capacity: 50ml · Joint: B-19, B-19, B-24 (centre, side necks)	1 pc	400
GLT-GW01-028	R.B Flask Three Neck With Joint — Capacity: 100ml · Joint: B-24, B-24, B-24 (centre, side necks)	1 pc	475
GLT-GW01-029	R.B Flask Three Neck With Joint — Capacity: 250ml · Joint: B-24, B-24, B-24 (centre, side necks)	1 pc	700
GLT-GW01-030	R.B Flask Three Neck With Joint — Capacity: 500ml · Joint: B-24, B-24, B-24 (centre, side necks)	1 pc	975
GLT-GW01-031	R.B Flask Three Neck With Joint — Capacity: 1000ml · Joint: B-24, B-24, B-24 (centre, side necks)	1 pc	1,275
GLT-GW01-032	R.B Flask Three Neck With Joint — Capacity: 2000ml · Joint: B-34, B-24, B-24 (centre, side necks)	1 pc	2,100
GLT-GW01-033	R.B Flask Three Neck With Joint — Capacity: 3000ml · Joint: B-34, B-24, B-24 (centre, side necks)	1 pc	3,525
GLT-GW01-034	R.B Flask Three Neck With Joint — Capacity: 5000ml · Joint: B-34, B-24, B-24 (centre, side necks)	1 pc	5,400
GLT-GW01-035	R.B Flask Three Neck With Joint — Capacity: 10000ml · Joint: B-45, B-24, B-24 (centre, side necks)	1 pc	10,250

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-04

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

R.B Flask Three Neck With Joint & R.B Flask Four Neck With Joint

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

7
PRODUCTS

7
PACK
OPTIONS

₹325
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
R.B Flask Three Neck With Joint			
GLT-GW01-036	R.B Flask Three Neck With Joint — Capacity: 20000ml · Joint: B-55, B-24, B-24 (centre, side necks)	1 pc	15,250
R.B Flask Four Neck With Joint			
GLT-GW01-037	R.B Flask Four Neck With Joint — Capacity: 25ml · Joint: B-14, B-14, B-14, B-14 (centre, side necks)	1 pc	325
GLT-GW01-038	R.B Flask Four Neck With Joint — Capacity: 50ml · Joint: B-19, B-19, B-14, B-14 (centre, side necks)	1 pc	475
GLT-GW01-039	R.B Flask Four Neck With Joint — Capacity: 100ml · Joint: B-24, B-24, B-24, B-24 (centre, side necks)	1 pc	550
GLT-GW01-040	R.B Flask Four Neck With Joint — Capacity: 250ml · Joint: B-24, B-24, B-24, B-24 (centre, side necks)	1 pc	775
GLT-GW01-041	R.B Flask Four Neck With Joint — Capacity: 500ml · Joint: B-24, B-24, B-24, B-24 (centre, side necks)	1 pc	1,075
GLT-GW01-042	R.B Flask Four Neck With Joint — Capacity: 1000ml · Joint: B-24, B-24, B-24, B-24 (centre, side necks)	1 pc	1,525

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-05

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

R.B Flask Four Neck With Joint & Beakers, low form, graduated with spout, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹109
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
R.B Flask Four Neck With Joint			
GLT-GW01-043	R.B Flask Four Neck With Joint — Capacity: 2000ml · Joint: B-34, B-24, B-24, B-24 (centre, side necks)	1 pc	2,325
GLT-GW01-044	R.B Flask Four Neck With Joint — Capacity: 3000ml · Joint: B-34, B-24, B-24, B-24 (centre, side necks)	1 pc	4,000
GLT-GW01-045	R.B Flask Four Neck With Joint — Capacity: 5000ml · Joint: B-34, B-24, B-24, B-24 (centre, side necks)	1 pc	5,750
GLT-GW01-046	R.B Flask Four Neck With Joint — Capacity: 10000ml · Joint: B-45, B-24, B-24, B-24 (centre, side necks)	1 pc	11,000
GLT-GW01-047	R.B Flask Four Neck With Joint — Capacity: 20000ml · Joint: B-55, B-24, B-24, B-24 (centre, side necks)	1 pc	18,000
Beakers, low form, graduated with spout, 3.3 borosilicate glass			
GLT-GW01-048	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 5* · O.D (mm): 20 · Height (mm): 25	1 pc (case of 20)	109

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-06

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Beakers, low form, graduated with spout, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

7
PRODUCTS

7
PACK
OPTIONS

₹109
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Beakers, low form, graduated with spout, 3.3 borosilicate glass			
GLT-GW01-049	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 10* · O.D (mm): 25 · Height (mm): 32	1 pc (case of 20)	112
GLT-GW01-050	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 25 · O.D (mm): 35 · Height (mm): 45	1 pc (case of 10)	125
GLT-GW01-051	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 50 · O.D (mm): 42 · Height (mm): 55	1 pc (case of 10)	109
GLT-GW01-052	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 100 · O.D (mm): 50 · Height (mm): 75	1 pc (case of 10)	119
GLT-GW01-053	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 150 · O.D (mm): 55 · Height (mm): 83	1 pc (case of 12)	125
GLT-GW01-054	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 250 · O.D (mm): 69 · Height (mm): 94	1 pc (case of 6)	125
GLT-GW01-055	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 500 · O.D (mm): 87 · Height (mm): 115	1 pc (case of 6)	200

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-07

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Beakers, low form, graduated with spout, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹204
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Beakers, low form, graduated with spout, 3.3 borosilicate glass			
GLT-GW01-056	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 600 · O.D (mm): 90 · Height (mm): 120	1 pc (case of 6)	204
GLT-GW01-057	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 1000 · O.D (mm): 110 · Height (mm): 150	1 pc (case of 2)	358
GLT-GW01-058	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 2000 · O.D (mm): 131 · Height (mm): 190	1 pc (case of 2)	684
GLT-GW01-059	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 3000 · O.D (mm): 146 · Height (mm): 214	1 pc	1,698
GLT-GW01-060	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 5000 · O.D (mm): 170 · Height (mm): 265	1 pc	2,462
GLT-GW01-061	Beakers, low form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 10000 · O.D (mm): 220 · Height (mm): 320	1 pc	8,078

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-08

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Beakers, tall form, graduated with spout, 3.3 borosilicate glass & Flasks, boiling, flat bottom, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹138
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Beakers, tall form, graduated with spout, 3.3 borosilicate glass			
GLT-GW01-062	Beakers, tall form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 50 · O.D (mm): 38 · Height (mm): 70	1 pc (case of 20)	178
GLT-GW01-063	Beakers, tall form, graduated with spout, 3.3 borosilicate glass — Capacity (ml): 100 · O.D (mm): 48 · Height (mm): 80	1 pc (case of 20)	181
Flasks, boiling, flat bottom, 3.3 borosilicate glass			
GLT-GW01-064	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 25 · Height (mm): 41 × 90 · O.D (mm): 20	1 pc (case of 2)	138
GLT-GW01-065	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 50 · Height (mm): 52 × 90 · O.D (mm): 26	1 pc (case of 2)	162
GLT-GW01-066	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 100 · Height (mm): 63 × 105 · O.D (mm): 26	1 pc (case of 2)	175
GLT-GW01-067	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 150 · Height (mm): 75 × 132 · O.D (mm): 28	1 pc (case of 2)	188

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-09

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, boiling, flat bottom, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

7
PRODUCTS

7
PACK
OPTIONS

₹206
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, boiling, flat bottom, 3.3 borosilicate glass			
GLT-GW01-068	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 250 · Height (mm): 85 × 140 · O.D (mm): 34	1 pc (case of 2)	206
GLT-GW01-069	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 300 · Height (mm): 87 × 160 · O.D (mm): 24	1 pc (case of 2)	256
GLT-GW01-070	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 500 · Height (mm): 105 × 170 · O.D (mm): 33	1 pc (case of 2)	269
GLT-GW01-071	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 1000 · Height (mm): 135 × 200 · O.D (mm): 42	1 pc (case of 2)	388
GLT-GW01-072	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 2000 · Height (mm): 167 × 250 · O.D (mm): 50	1 pc	712
GLT-GW01-073	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 3000 · Height (mm): 186 × 300 · O.D (mm): 50	1 pc	2,100
GLT-GW01-074	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 5000 · Height (mm): 223 × 340 · O.D (mm): 60	1 pc	3,200

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-10

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, boiling, flat bottom, 3.3 borosilicate glass & Flasks, boiling, round bottom, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹138
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, boiling, flat bottom, 3.3 borosilicate glass			
GLT-GW01-075	Flasks, boiling, flat bottom, 3.3 borosilicate glass — Capacity (ml): 10000 · Height (mm): 279 × 400 · O.D (mm): 65	1 pc	6,738
Flasks, boiling, round bottom, 3.3 borosilicate glass			
GLT-GW01-076	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 25 · Dimensions (mm): 41 × 90 · Approx. neck OD. (mm): 20	1 pc (case of 2)	138
GLT-GW01-077	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 50 · Dimensions (mm): 52 × 90 · Approx. neck OD. (mm): 26	1 pc (case of 2)	162
GLT-GW01-078	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 100 · Dimensions (mm): 63 × 105 · Approx. neck OD. (mm): 26	1 pc (case of 2)	194
GLT-GW01-079	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 150 · Dimensions (mm): 75 × 132 · Approx. neck OD. (mm): 28	1 pc (case of 2)	206
GLT-GW01-080	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 250 · Dimensions (mm): 85 × 140 · Approx. neck OD. (mm): 34	1 pc (case of 2)	256

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-11

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, boiling, round bottom, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹269
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, boiling, round bottom, 3.3 borosilicate glass			
GLT-GW01-081	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 500 · Dimensions (mm): 105 × 170 · Approx. neck OD. (mm): 33	1 pc (case of 2)	269
GLT-GW01-082	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 1000 · Dimensions (mm): 135 × 200 · Approx. neck OD. (mm): 42	1 pc (case of 2)	388
GLT-GW01-083	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 2000 · Dimensions (mm): 167 × 250 · Approx. neck OD. (mm): 50	1 pc (case of 2)	712
GLT-GW01-084	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 3000 · Dimensions (mm): 186 × 300 · Approx. neck OD. (mm): 50	1 pc	2,100
GLT-GW01-085	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 5000 · Dimensions (mm): 223 × 340 · Approx. neck OD. (mm): 60	1 pc	3,200
GLT-GW01-086	Flasks, boiling, round bottom, 3.3 borosilicate glass — Capacity (ml): 10000 · Dimensions (mm): 279 × 400 · Approx. neck OD. (mm): 65	1 pc	6,738

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-12

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, boiling, flat bottom, short neck with I/C joint, 3.3 borosilicate glass & Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹238
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, boiling, flat bottom, short neck with I/C joint, 3.3 borosilicate glass			
GLT-GW01-087	Flasks, boiling, flat bottom, short neck with I/C joint, 3.3 borosilicate glass — Capacity (ml): 150 · Joint: 24/29 · Height (mm): 115	1 pc (case of 2)	256
GLT-GW01-088	Flasks, boiling, flat bottom, short neck with I/C joint, 3.3 borosilicate glass — Capacity (ml): 250 · Joint: 24/29 · Height (mm): 130	1 pc (case of 2)	300
GLT-GW01-089	Flasks, boiling, flat bottom, short neck with I/C joint, 3.3 borosilicate glass — Capacity (ml): 500 · Joint: 24/29 · Height (mm): 148	1 pc (case of 2)	400
Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass			
GLT-GW01-090	Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass — Capacity (ml): 50 · Joint: 19/26 · Height (mm): 90	1 pc (case of 2)	238
GLT-GW01-091	Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass — Capacity (ml): 100 · Joint: 24/29 · Height (mm): 110	1 pc (case of 2)	256
GLT-GW01-092	Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass — Capacity (ml): 150 · Joint: 24/29 · Height (mm): 115	1 pc (case of 2)	288

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-13

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass & Flasks, boiling, round bottom, two necks, centre and one angled side neck, I/C joint, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹319
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass			
GLT-GW01-093	Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass — Capacity (ml): 250 · Joint: 24/29 · Height (mm): 130	1 pc (case of 2)	319
GLT-GW01-094	Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass — Capacity (ml): 500 · Joint: 24/29 · Height (mm): 150	1 pc (case of 2)	400
GLT-GW01-095	Flasks, boiling, round bottom, short neck, with I/C joint, 3.3 borosilicate glass — Capacity (ml): 1000 · Joint: 24/29 · Height (mm): 188	1 pc (case of 2)	550
Flasks, boiling, round bottom, two necks, centre and one angled side neck, I/C joint, 3.3 borosilicate glass			
GLT-GW01-096	Flasks, boiling, round bottom, two necks, centre and one angled side neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 100 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	375
GLT-GW01-097	Flasks, boiling, round bottom, two necks, centre and one angled side neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 250 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	881
GLT-GW01-098	Flasks, boiling, round bottom, two necks, centre and one angled side neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 500 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	588

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-14

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, boiling, round bottom, two necks, centre and one angled side neck, I/C joint, 3.3 borosilicate glass & Flasks, boiling, round bottom, two necks, centre and one parallel side neck, I/C joint, 3.3 borosilicate glass

Illustrative image

Low - and tall - form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

5
PRODUCTS

5
PACK
OPTIONS

₹331
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, boiling, round bottom, two necks, centre and one angled side neck, I/C joint, 3.3 borosilicate glass			
GLT-GW01-099	Flasks, boiling, round bottom, two necks, centre and one angled side neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 1000 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	838
Flasks, boiling, round bottom, two necks, centre and one parallel side neck, I/C joint, 3.3 borosilicate glass			
GLT-GW01-100	Flasks, boiling, round bottom, two necks, centre and one parallel side neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 100 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	331
GLT-GW01-101	Flasks, boiling, round bottom, two necks, centre and one parallel side neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 250 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	469
GLT-GW01-102	Flasks, boiling, round bottom, two necks, centre and one parallel side neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 500 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	650
GLT-GW01-103	Flasks, boiling, round bottom, two necks, centre and one parallel side neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 1000 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	812

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-15

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, boiling, round bottom, three necks, centre and two angled side necks, I/C joint, 3.3 borosilicate glass & Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

5
PRODUCTS

5
PACK
OPTIONS

₹438
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, boiling, round bottom, three necks, centre and two angled side necks, I/C joint, 3.3 borosilicate glass			
GLT-GW01-104	Flasks, boiling, round bottom, three necks, centre and two angled side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 100 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	438
GLT-GW01-105	Flasks, boiling, round bottom, three necks, centre and two angled side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 250 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	606
GLT-GW01-106	Flasks, boiling, round bottom, three necks, centre and two angled side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 500 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	881
GLT-GW01-107	Flasks, boiling, round bottom, three necks, centre and two angled side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 1000 · Joint: 24/29 · Joint: 19/26	1 pc (case of 10)	1,112
Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass			
GLT-GW01-108	Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 100 · Joint: 24/29 · Joint: 19/26	1 pc	438

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-16

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

5
PRODUCTS

5
PACK
OPTIONS

₹606
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass			
GLT-GW01-109	Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 250 · Joint: 24/29 · Joint: 19/26	1 pc	606
GLT-GW01-110	Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 500 · Joint: 24/29 · Joint: 19/26	1 pc	881
GLT-GW01-111	Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 1000 · Joint: 24/29 · Joint: 19/26	1 pc	1,112
GLT-GW01-112	Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 2000 · Joint: 24/29 · Joint: 19/26	1 pc	2,300
GLT-GW01-113	Flasks, boiling, round bottom, three necks, centre and two parallel side necks, I/C joint, 3.3 borosilicate glass — Capacity (ml): 5000 · Joint: 24/29 · Joint: 19/26	1 pc	7,525

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-17

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Erlenmeyer, conical, with I/C joint, 3.3 borosilicate glass & Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹300
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Erlenmeyer, conical, with I/C joint, 3.3 borosilicate glass			
GLT-GW01-114	Flasks, Erlenmeyer, conical, with I/C joint, 3.3 borosilicate glass — Capacity (ml): 100 · Height (mm): 90 · Joint: 24/25	1 pc (case of 2)	300
GLT-GW01-115	Flasks, Erlenmeyer, conical, with I/C joint, 3.3 borosilicate glass — Capacity (ml): 250 · Height (mm): 145 · Joint: 24/29	1 pc (case of 2)	362
GLT-GW01-116	Flasks, Erlenmeyer, conical, with I/C joint, 3.3 borosilicate glass — Capacity (ml): 500 · Height (mm): 175 · Joint: 24/29	1 pc (case of 2)	500
Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass			
GLT-GW01-117	Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 25 · Height (mm): 65 · Joint: 14/15	1 pc (case of 2)	312
GLT-GW01-118	Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 50 · Height (mm): 85 · Joint: 19/20	1 pc (case of 2)	344
GLT-GW01-119	Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 100 · Height (mm): 90 · Joint: 24/25	1 pc (case of 2)	388

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-18

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass & Flask, Conical, with Screw Cap, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

8
PRODUCTS

8
PACK
OPTIONS

₹300
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass			
GLT-GW01-120	Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 250 · Height (mm): 145 · Joint: 29/25	1 pc (case of 2)	438
GLT-GW01-121	Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 500 · Height (mm): 175 · Joint: 29/25	1 pc (case of 2)	575
GLT-GW01-122	Flasks, Erlenmeyer, conical, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 1000 · Height (mm): 175 · Joint: 29/25	1 pc (case of 2)	962
Flask, Conical, with Screw Cap, 3.3 borosilicate glass			
GLT-GW01-123	Flask, Conical, with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 100 · Height (mm): 100	1 pc (case of 20)	300
GLT-GW01-124	Flask, Conical, with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 150 · Height (mm): 125	1 pc (case of 20)	362
GLT-GW01-125	Flask, Conical, with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 250 · Height (mm): 145	1 pc (case of 20)	450
GLT-GW01-126	Flask, Conical, with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 500 · Height (mm): 180	1 pc (case of 10)	550
GLT-GW01-127	Flask, Conical, with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 1000 · Height (mm): 225	1 pc (case of 10)	881

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-19

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flask, Conical, Amber with Screw Cap, 3.3 borosilicate glass & Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

7
PRODUCTS

7
PACK
OPTIONS

₹138
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flask, Conical, Amber with Screw Cap, 3.3 borosilicate glass			
GLT-GW01-128	Flask, Conical, Amber with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 100 · Height (mm): 100	1 pc (case of 20)	619
GLT-GW01-129	Flask, Conical, Amber with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 150 · Height (mm): 125	1 pc (case of 20)	775
GLT-GW01-130	Flask, Conical, Amber with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 250 · Height (mm): 145	1 pc (case of 20)	881
GLT-GW01-131	Flask, Conical, Amber with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 500 · Height (mm): 180	1 pc (case of 10)	1,088
GLT-GW01-132	Flask, Conical, Amber with Screw Cap, 3.3 borosilicate glass — Capacity (ml): 1000 · Height (mm): 225	1 pc (case of 10)	2,488
Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass			
GLT-GW01-133	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 10 · Dimensions (mm): 30 × 55 · Neck O.D. (mm): 15	1 pc (case of 12)	138
GLT-GW01-134	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 25 · Dimensions (mm): 40 × 70 · Neck O.D. (mm): 16	1 pc (case of 12)	138

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-20

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

7
PRODUCTS

7
PACK
OPTIONS

₹138
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass			
GLT-GW01-135	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 50 · Dimensions (mm): 52 x 90 · Neck O.D. (mm): 20	1 pc (case of 12)	138
GLT-GW01-136	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 100 · Dimensions (mm): 60 x 105 · Neck O.D. (mm): 23	1 pc (case of 2)	138
GLT-GW01-137	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 150 · Dimensions (mm): 69 x 120 · Neck O.D. (mm): 25	1 pc (case of 2)	188
GLT-GW01-138	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 250 · Dimensions (mm): 83 x 184 · Neck O.D. (mm): 30	1 pc (case of 6)	225
GLT-GW01-139	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 500 · Dimensions (mm): 100 x 195 · Neck O.D. (mm): 35	1 pc (case of 6)	300
GLT-GW01-140	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 1000 · Dimensions (mm): 128 x 225 · Neck O.D. (mm): 40	1 pc (case of 6)	425
GLT-GW01-141	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 2000 · Dimensions (mm): 175 x 285 · Neck O.D. (mm): 48	1 pc (case of 6)	775

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-21

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass & Flasks, Erlenmeyer, graduated, conical, Amber, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹269
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass			
GLT-GW01-142	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 3000 · Dimensions (mm): 185 x 315 · Neck O.D. (mm): 48	1 pc (case of 2)	2,100
GLT-GW01-143	Flasks, Erlenmeyer, graduated, conical, 3.3 borosilicate glass — Capacity (ml): 5000 · Dimensions (mm): 220 x 377 · Neck O.D. (mm): 50	1 pc (case of 4)	3,200
Flasks, Erlenmeyer, graduated, conical, Amber, 3.3 borosilicate glass			
GLT-GW01-144	Flasks, Erlenmeyer, graduated, conical, Amber, 3.3 borosilicate glass — Capacity (ml): 100 · Dimensions (mm): 60 × 105 · Neck O.D. (mm): 23	1 pc (case of 2)	269
GLT-GW01-145	Flasks, Erlenmeyer, graduated, conical, Amber, 3.3 borosilicate glass — Capacity (ml): 150 · Dimensions (mm): 69 × 120 · Neck O.D. (mm): 25	1 pc (case of 2)	312
GLT-GW01-146	Flasks, Erlenmeyer, graduated, conical, Amber, 3.3 borosilicate glass — Capacity (ml): 250 · Dimensions (mm): 83 × 184 · Neck O.D. (mm): 30	1 pc (case of 2)	394
GLT-GW01-147	Flasks, Erlenmeyer, graduated, conical, Amber, 3.3 borosilicate glass — Capacity (ml): 500 · Dimensions (mm): 100 × 195 · Neck O.D. (mm): 35	1 pc (case of 2)	525

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-22

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Erlenmeyer, graduated, conical, Amber, 3.3 borosilicate glass, Flasks, Filtration, Heavy Wall, with outlet for Vacuum, 3.3 borosilicate glass & more

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹425
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Erlenmeyer, graduated, conical, Amber, 3.3 borosilicate glass			
GLT-GW01-148	Flasks, Erlenmeyer, graduated, conical, Amber, 3.3 borosilicate glass — Capacity (ml): 1000 · Dimensions (mm): 128 × 225 · Neck O.D. (mm): 40	1 pc (case of 2)	950
Flasks, Filtration, Heavy Wall, with outlet for Vacuum, 3.3 borosilicate glass			
GLT-GW01-149	Flasks, Filtration, Heavy Wall, with outlet for Vacuum, 3.3 borosilicate glass — Capacity (ml): 500 · Dimensions (mm): 105 × 184	1 pc (case of 2)	725
GLT-GW01-150	Flasks, Filtration, Heavy Wall, with outlet for Vacuum, 3.3 borosilicate glass — Capacity (ml): 1000 · Dimensions (mm): 136 × 250	1 pc (case of 2)	1,075
GLT-GW01-151	Flasks, Filtration, Heavy Wall, with outlet for Vacuum, 3.3 borosilicate glass — Capacity (ml): 2000 · Dimensions (mm): 168 × 295	1 pc (case of 2)	2,200
GLT-GW01-152	Flasks, Filtration, Heavy Wall, with outlet for Vacuum, 3.3 borosilicate glass — Capacity (ml): 3000 · Dimensions (mm): 237 × 450	1 pc	8,044
Flasks, Iodine determination, with stopper, 3.3 borosilicate glass			
GLT-GW01-153	Flasks, Iodine determination, with stopper, 3.3 borosilicate glass — Capacity (ml): 250 · Joint: 29/25	1 pc (case of 2)	425

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-23

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Iodine determination, with stopper, 3.3 borosilicate glass, Flasks, Kjeldahl, round bottom, long neck, 3.3 borosilicate glass & more

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹538
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Iodine determination, with stopper, 3.3 borosilicate glass			
GLT-GW01-154	Flasks, Iodine determination, with stopper, 3.3 borosilicate glass — Capacity (ml): 500 · Joint: 29/25	1 pc (case of 2)	538
Flasks, Kjeldahl, round bottom, long neck, 3.3 borosilicate glass			
GLT-GW01-155	Flasks, Kjeldahl, round bottom, long neck, 3.3 borosilicate glass — Capacity (ml): 300 · O.D (mm): 80 · Height (mm): 300	1 pc (case of 2)	538
GLT-GW01-156	Flasks, Kjeldahl, round bottom, long neck, 3.3 borosilicate glass — Capacity (ml): 500 · O.D (mm): 97 · Height (mm): 330	1 pc (case of 2)	850
GLT-GW01-157	Flasks, Kjeldahl, round bottom, long neck, 3.3 borosilicate glass — Capacity (ml): 800 · O.D (mm): 114 · Height (mm): 378	1 pc	1,156
Flasks, Kjeldahl, round bottom, long neck, I/C joint, 3.3 borosilicate glass			
GLT-GW01-158	Flasks, Kjeldahl, round bottom, long neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 500 · O.D (mm): 97 · Height (mm): 330	1 pc (case of 2)	1,031
GLT-GW01-159	Flasks, Kjeldahl, round bottom, long neck, I/C joint, 3.3 borosilicate glass — Capacity (ml): 800 · O.D (mm): 114 · Height (mm): 378	1 pc	1,112

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-24

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Volumetric, Sugar Estimation, without stopper, 3.3 borosilicate glass & Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate

Illustrative image

Low - and tall - form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

5
PRODUCTS

5
PACK
OPTIONS

₹512
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Volumetric, Sugar Estimation, without stopper, 3.3 borosilicate glass			
GLT-GW01-160	Flasks, Volumetric, Sugar Estimation, without stopper, 3.3 borosilicate glass — aan: 100/110 · 0.15	1 pc (case of 2)	512
Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate			
GLT-GW01-161	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 5 · Tolerance ($\pm$ ml): 0.02 · Stopper size: 10//15	1 pc (case of 2)	1,125
GLT-GW01-162	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 10 · Tolerance ($\pm$ ml): 0.02 · Stopper size: 10//15	1 pc (case of 2)	975
GLT-GW01-163	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 20 · Tolerance ($\pm$ ml): 0.03 · Stopper size: 10//15	1 pc (case of 2)	975
GLT-GW01-164	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 25 · Tolerance ($\pm$ ml): 0.03 · Stopper size: 10//15	1 pc (case of 2)	944

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-25

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate

Illustrative image

Low - and tall - form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

5
PRODUCTS

5
PACK
OPTIONS

₹925
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate			
GLT-GW01-165	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 50 · Tolerance ($\pm$ ml): 0.05 · Stopper size: 10//15	1 pc (case of 2)	944
GLT-GW01-166	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 100 · Tolerance ($\pm$ ml): 0.08 · Stopper size: 14//15	1 pc (case of 2)	925
GLT-GW01-167	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 200 · Tolerance ($\pm$ ml): 0.1 · Stopper size: 14//15	1 pc (case of 2)	1,262
GLT-GW01-168	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 250 · Tolerance ($\pm$ ml): 0.12 · Stopper size: 14//15	1 pc (case of 2)	1,262
GLT-GW01-169	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 500 · Tolerance ($\pm$ ml): 0.15 · Stopper size: 19//20	1 pc (case of 2)	1,675

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-26

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate & Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass

Illustrative image

Low - and tall - form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹538
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate			
GLT-GW01-170	Flask, Volumetric, with stopper, 3.3 borosilicate glass, As per USP with individual certificate — Capacity (ml): 1000 · Tolerance (± ml): 0.3 · Stopper size: 19//20	1 pc	2,406
Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass			
GLT-GW01-171	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 5 · Tolerance (± ml): 0.05 · Stopper size: 10/15	1 pc (case of 2)	606
GLT-GW01-172	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 10 · Tolerance (± ml): 0.05 · Stopper size: 10/15	1 pc (case of 2)	606
GLT-GW01-173	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 20 · Tolerance (± ml): 0.05 · Stopper size: 10/15	1 pc (case of 2)	538
GLT-GW01-174	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 25 · Tolerance (± ml): 0.08 · Stopper size: 10/15	1 pc (case of 2)	575
GLT-GW01-175	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 50 · Tolerance (± ml): 0.12 · Stopper size: 10/15	1 pc (case of 2)	588

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-27

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹606
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass			
GLT-GW01-176	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 100 · Tolerance ($\pm$ ml): 0.20 · Stopper size: 14/15	1 pc (case of 2)	606
GLT-GW01-177	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 200 · Tolerance ($\pm$ ml): 0.30 · Stopper size: 14/15	1 pc (case of 2)	744
GLT-GW01-178	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 250 · Tolerance ($\pm$ ml): 0.30 · Stopper size: 14/15	1 pc (case of 2)	625
GLT-GW01-179	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 500 · Tolerance ($\pm$ ml): 0.50 · Stopper size: 19/20	1 pc (case of 2)	962
GLT-GW01-180	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 1000 · Tolerance ($\pm$ ml): 0.80 · Stopper size: 19/20	1 pc	1,612
GLT-GW01-181	Flasks, Volumetric, Amber, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 2000 · Tolerance ($\pm$ ml): 0.80 · Stopper size: 19/20	1 pc	2,600

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-28

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A'

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

5
PRODUCTS

5
PACK
OPTIONS

₹1,062
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate)			
GLT-GW01-182	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 5 · Tolerance ($\pm$ ml): 0.025 · Stopper size: 10/15	1 pc (case of 2)	1,250
GLT-GW01-183	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 10 · Tolerance ($\pm$ ml): 0.025 · Stopper size: 10/15	1 pc (case of 2)	1,169
GLT-GW01-184	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 20 · Tolerance ($\pm$ ml): 0.04 · Stopper size: 10/15	1 pc (case of 2)	1,169
GLT-GW01-185	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 25 · Tolerance ($\pm$ ml): 0.04 · Stopper size: 10/15	1 pc (case of 2)	1,088
GLT-GW01-186	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 50 · Tolerance ($\pm$ ml): 0.06 · Stopper size: 10/15	1 pc (case of 2)	1,062

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-29

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A'

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

5
PRODUCTS

5
PACK
OPTIONS

₹1,112
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate)			
GLT-GW01-187	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 100 · Tolerance ($\pm$ ml): 0.10 · Stopper size: 14/15	1 pc (case of 2)	1,112
GLT-GW01-188	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 200 · Tolerance ($\pm$ ml): 0.15 · Stopper size: 14/15	1 pc (case of 2)	1,562
GLT-GW01-189	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 250 · Tolerance ($\pm$ ml): 0.15 · Stopper size: 14/15	1 pc (case of 2)	1,569
GLT-GW01-190	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 500 · Tolerance ($\pm$ ml): 0.25 · Stopper size: 19/20	1 pc (case of 2)	2,031
GLT-GW01-191	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 1000 · Tolerance ($\pm$ ml): 0.40 · Stopper size: 19/20	1 pc	2,988

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-30

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' & Flasks, Volumetric, I/C stopper, Class 'A', 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

5
PRODUCTS

5
PACK
OPTIONS

₹850
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate)			
GLT-GW01-192	Flasks, Volumetric, with I/C stopper, Amber, 3.3 borosilicate glass, Class 'A' (with certificate) — Capacity (ml): 2000 · Tolerance ($\pm$ ml): 0.60 · Stopper size: 24/25	1 pc	4,044
Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass			
GLT-GW01-193	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 5 · Tolerance ($\pm$ ml): 0.025 · Stopper size: 10/15	1 pc (case of 2)	1,125
GLT-GW01-194	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 10 · Tolerance ($\pm$ ml): 0.025 · Stopper size: 10/15	1 pc (case of 2)	950
GLT-GW01-195	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 20 · Tolerance ($\pm$ ml): 0.04 · Stopper size: 10/15	1 pc (case of 2)	950
GLT-GW01-196	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 25 · Tolerance ($\pm$ ml): 0.04 · Stopper size: 10/15	1 pc (case of 2)	850

Full specifications and every pack size: Product Catalogue 2026-27, pages **470-480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW01-31

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Volumetric, I/C stopper, Class 'A', 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

5
PRODUCTS

5
PACK
OPTIONS

₹838
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass			
GLT-GW01-197	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 50 · Tolerance ($\pm$ ml): 0.06 · Stopper size: 10/15	1 pc (case of 2)	850
GLT-GW01-198	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 100 · Tolerance ($\pm$ ml): 0.10 · Stopper size: 14/15	1 pc (case of 2)	838
GLT-GW01-199	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 200 · Tolerance ($\pm$ ml): 0.15 · Stopper size: 14/15	1 pc (case of 2)	1,050
GLT-GW01-200	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 250 · Tolerance ($\pm$ ml): 0.15 · Stopper size: 14/15	1 pc (case of 2)	1,125
GLT-GW01-201	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 500 · Tolerance ($\pm$ ml): 0.25 · Stopper size: 19/20	1 pc (case of 2)	1,500

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-32

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Volumetric, I/C stopper, Class 'A', 3.3 borosilicate glass & Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

6
PRODUCTS

6
PACK
OPTIONS

₹331
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass			
GLT-GW01-202	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 1000 · Tolerance ($\pm$ ml): 0.40 · Stopper size: 19/20	1 pc	2,225
GLT-GW01-203	Flasks, Volumetric, I/C stopper, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 2000 · Tolerance ($\pm$ ml): 0.60 · Stopper size: 24/25	1 pc	3,438
Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass			
GLT-GW01-204	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 20 · Tolerance ($\pm$ ml): 0.05 · Stopper size: 10/15	1 pc (case of 2)	331
GLT-GW01-205	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 5 · Tolerance ($\pm$ ml): 0.05 · Stopper size: 10/15	1 pc (case of 2)	362
GLT-GW01-206	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 10 · Tolerance ($\pm$ ml): 0.05 · Stopper size: 10/15	1 pc (case of 2)	344
GLT-GW01-207	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 25 · Tolerance ($\pm$ ml): 0.08 · Stopper size: 10/15	1 pc (case of 2)	344

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-33

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW01 GLASS BEAKERS & FLASKS

Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass

Illustrative image

Low - and tall-form beakers and a wide choice of flasks — conical, volumetric, boiling, distilling, Kjeldahl and filtering — mostly in 3.3 borosilicate glass.

7
PRODUCTS

7
PACK
OPTIONS

₹362
PRICE FROM

TYPICAL RESEARCH USE

Solution preparation and mixing

Titration and volumetric work

Heating, boiling and distillation set-ups

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass			
GLT-GW01-208	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 50 · Tolerance ($\pm$ ml): 0.12 · Stopper size: 10/15	1 pc (case of 2)	362
GLT-GW01-209	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 100 · Tolerance ($\pm$ ml): 0.20 · Stopper size: 14/15	1 pc (case of 2)	425
GLT-GW01-210	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 200 · Tolerance ($\pm$ ml): 0.30 · Stopper size: 14/15	1 pc (case of 2)	475
GLT-GW01-211	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 250 · Tolerance ($\pm$ ml): 0.30 · Stopper size: 14/15	1 pc (case of 2)	500
GLT-GW01-212	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 500 · Tolerance ($\pm$ ml): 0.50 · Stopper size: 19/20	1 pc (case of 2)	638
GLT-GW01-213	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 1000 · Tolerance ($\pm$ ml): 0.80 · Stopper size: 19/20	1 pc	950
GLT-GW01-214	Flasks, Volumetric, I/C stopper, 3.3 borosilicate glass — Capacity (ml): 2000 · Tolerance ($\pm$ ml): 1.20 · Stopper size: 24/25	1 pc	1,800

Full specifications and every pack size: Product Catalogue 2026–27, pages **470–480**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW01-34

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113