

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Burettes, straight bore stopcock with Glass key, Class 'A', amber stain graduation & Burettes, straight bore stopcock with PTFE key, Class 'A', amber stain graduation

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹2,250
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Burettes, straight bore stopcock with Glass key, Class 'A' (with certificate), amber stain graduation			
GLT-GW03-001	Burettes, straight bore stopcock with Glass key, Class 'A' (with certificate), amber stain graduation — Capacity (ml): 10 · Graduation (ml): 0.05 · Tolerance ($\pm$ ml): 0.02	1 pc (case of 5)	2,250
GLT-GW03-002	Burettes, straight bore stopcock with Glass key, Class 'A' (with certificate), amber stain graduation — Capacity (ml): 25 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.05	1 pc (case of 5)	2,338
GLT-GW03-003	Burettes, straight bore stopcock with Glass key, Class 'A' (with certificate), amber stain graduation — Capacity (ml): 50 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.05	1 pc (case of 5)	2,675
GLT-GW03-004	Burettes, straight bore stopcock with Glass key, Class 'A' (with certificate), amber stain graduation — Capacity (ml): 100 · Graduation (ml): 0.20 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 5)	3,800
Burettes, straight bore stopcock with PTFE key, Class 'A' (with certificate), amber stain graduation			
GLT-GW03-005	Burettes, straight bore stopcock with PTFE key, Class 'A' (with certificate), amber stain graduation — Capacity (ml): 10 · Graduation (ml): 0.05 · Tolerance ($\pm$ ml): 0.02	1 pc (case of 5)	2,250

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-01

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Burettes, straight bore stopcock with PTFE key, Class 'A', amber stain graduation & Burettes, straight bore stopcock with Glass key, amber stain graduation, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹844
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Burettes, straight bore stopcock with PTFE key, Class 'A' (with certificate), amber stain graduation			
GLT-GW03-006	Burettes, straight bore stopcock with PTFE key, Class 'A' (with certificate), amber stain graduation — Capacity (ml): 25 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.05	1 pc (case of 5)	2,425
GLT-GW03-007	Burettes, straight bore stopcock with PTFE key, Class 'A' (with certificate), amber stain graduation — Capacity (ml): 50 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.05	1 pc (case of 5)	2,762
GLT-GW03-008	Burettes, straight bore stopcock with PTFE key, Class 'A' (with certificate), amber stain graduation — Capacity (ml): 100 · Graduation (ml): 0.20 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 5)	3,450
Burettes, straight bore stopcock with Glass key, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-009	Burettes, straight bore stopcock with Glass key, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.05 · Tolerance ($\pm$ ml): 0.05	1 pc (case of 12)	844
GLT-GW03-010	Burettes, straight bore stopcock with Glass key, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 25 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 12)	844

Full specifications and every pack size: Product Catalogue 2026-27, pages **485-491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW03-02

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Burettes, straight bore stopcock with Glass key, amber stain graduation, 3.3 borosilicate glass, Burettes, straight bore stopcock with PTFE key, white stain graduation, 3.3 borosilicate glass & more

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹750
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Burettes, straight bore stopcock with Glass key, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-011	Burettes, straight bore stopcock with Glass key, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 50 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 12)	912
GLT-GW03-012	Burettes, straight bore stopcock with Glass key, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 100 · Graduation (ml): 0.20 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 12)	1,250
Burettes, straight bore stopcock with PTFE key, white stain graduation, 3.3 borosilicate glass			
GLT-GW03-013	Burettes, straight bore stopcock with PTFE key, white stain graduation, 3.3 borosilicate glass — Capacity (ml): 50 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 12)	1,100
Burettes, screw - thread (PTFE key) screw thread stopcock with PTFE key, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-014	Burettes, screw - thread (PTFE key) screw thread stopcock with PTFE key, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.05 · Tolerance ($\pm$ ml): 0.05	1 pc (case of 12)	750
GLT-GW03-015	Burettes, screw - thread (PTFE key) screw thread stopcock with PTFE key, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 25 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 12)	750

Full specifications and every pack size: Product Catalogue 2026-27, pages **485-491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW03-03

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Burettes, screw-thread screw thread stopcock with PTFE key, amber stain graduation, 3.3 borosilicate glass, Cylinders, Nessler, colour comparison, flat bottom, graduated, 3.3 borosilicate glass & more

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹250
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Burettes, screw - thread (PTFE key) screw thread stopcock with PTFE key, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-016	Burettes, screw - thread (PTFE key) screw thread stopcock with PTFE key, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 50 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 12)	788
GLT-GW03-017	Burettes, screw - thread (PTFE key) screw thread stopcock with PTFE key, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 100 · Graduation (ml): 0.20 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 12)	1,238
Cylinders, Nessler, colour comparison, flat bottom, graduated, 3.3 borosilicate glass			
GLT-GW03-018	Cylinders, Nessler, colour comparison, flat bottom, graduated, 3.3 borosilicate glass — Capacity (ml): 50 · O.D (mm): 26 · Height (mm): 150	1 pc (case of 80)	250
GLT-GW03-019	Cylinders, Nessler, colour comparison, flat bottom, graduated, 3.3 borosilicate glass — Capacity (ml): 100 · O.D (mm): 34 · Height (mm): 180	1 pc (case of 80)	350
Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-020	Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.2 · Tolerance ($\pm$ ml): 0.2 · Height (mm): 140	1 pc (case of 2)	281

Full specifications and every pack size: Product Catalogue 2026-27, pages **485-491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW03-04

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹300
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-021	Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 25 · Graduation (ml): 0.5 · Tolerance ($\pm$ ml): 0.5 · Height (mm): 170	1 pc (case of 2)	300
GLT-GW03-022	Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 50 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 1.0 · Height (mm): 200	1 pc (case of 2)	400
GLT-GW03-023	Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 100 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 1.0 · Height (mm): 260	1 pc (case of 2)	500
GLT-GW03-024	Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 250 · Graduation (ml): 2.0 · Tolerance ($\pm$ ml): 2.0 · Height (mm): 335	1 pc (case of 2)	750
GLT-GW03-025	Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 500 · Graduation (ml): 5.0 · Tolerance ($\pm$ ml): 5.0 · Height (mm): 390	1 pc	788

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-05

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass & Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A', 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

4
PRODUCTS

4
PACK
OPTIONS

₹688
PRICE FROM

TYPICAL RESEARCH USE

Titrations

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-026	Cylinders, graduated, with spout, plastic base and safety guard, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 1000 · Graduation (ml): 10.0 · Tolerance ($\pm$ ml): 10.0 · Height (mm): 470	1 pc	1,475
Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A' (with certificate), 3.3 borosilicate glass			
GLT-GW03-027	Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.2 · Tolerance ($\pm$ ml): 0.10 · Height (mm): 140	1 pc (case of 2)	688
GLT-GW03-028	Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 25 · Graduation (ml): 0.5 · Tolerance ($\pm$ ml): 0.25 · Height (mm): 170	1 pc (case of 2)	756
GLT-GW03-029	Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 50 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 0.50 · Height (mm): 200	1 pc (case of 2)	925

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-06

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A', 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

4
PRODUCTS

4
PACK
OPTIONS

₹1,031
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A' (with certificate), 3.3 borosilicate glass			
GLT-GW03-030	Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 100 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 0.50 · Height (mm): 260	1 pc (case of 2)	1,031
GLT-GW03-031	Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 250 · Graduation (ml): 2.0 · Tolerance ($\pm$ ml): 1.00 · Height (mm): 335	1 pc (case of 2)	1,556
GLT-GW03-032	Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 500 · Graduation (ml): 5.0 · Tolerance ($\pm$ ml): 2.50 · Height (mm): 390	1 pc (case of 2)	2,131
GLT-GW03-033	Cylinders, graduated, hexagonal base and spout, amber stain graduation, Class 'A' (with certificate), 3.3 borosilicate glass — Capacity (ml): 1000 · Graduation (ml): 10.0 · Tolerance ($\pm$ ml): 5.00 · Height (mm): 470	1 pc (case of 2)	2,544

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-07

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹375
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-034	Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.2 · Tolerance ($\pm$ ml): 0.2 · Height (mm): 140	1 pc (case of 2)	375
GLT-GW03-035	Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 25 · Graduation (ml): 0.5 · Tolerance ($\pm$ ml): 0.5 · Height (mm): 170	1 pc (case of 2)	425
GLT-GW03-036	Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 50 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 1.0 · Height (mm): 200	1 pc (case of 2)	481
GLT-GW03-037	Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 100 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 1.0 · Height (mm): 260	1 pc (case of 2)	575
GLT-GW03-038	Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 250 · Graduation (ml): 2.0 · Tolerance ($\pm$ ml): 2.0 · Height (mm): 335	1 pc (case of 2)	1,031

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-08

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass & Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹138
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-039	Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 500 · Graduation (ml): 5.0 · Tolerance ($\pm$ ml): 5.0 · Height (mm): 390	1 pc	138
GLT-GW03-040	Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 1000 · Graduation (ml): 10.0 · Tolerance ($\pm$ ml): 10.0 · Height (mm): 470	1 pc	2,025
GLT-GW03-041	Cylinders, graduated with round base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 2000 · Graduation (ml): 10.0 · Tolerance ($\pm$ ml): 10.0 · Height (mm): 570	1 pc	3,438
Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-042	Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.2 · Tolerance ($\pm$ ml): 0.2	1 pc (case of 2)	425
GLT-GW03-043	Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 25 · Graduation (ml): 0.5 · Tolerance ($\pm$ ml): 0.5	1 pc (case of 2)	475

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-09

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹500
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass			
GLT-GW03-044	Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 50 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 1.0	1 pc (case of 2)	500
GLT-GW03-045	Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 100 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 1.0	1 pc (case of 2)	575
GLT-GW03-046	Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 250 · Graduation (ml): 2.0 · Tolerance ($\pm$ ml): 2.0	1 pc	1,000
GLT-GW03-047	Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 500 · Graduation (ml): 5.0 · Tolerance ($\pm$ ml): 5.0	1 pc	1,375
GLT-GW03-048	Cylinders, measuring, with hexagonal base and spout, amber stain graduation, 3.3 borosilicate glass — Capacity (ml): 1000 · Graduation (ml): 10.0 · Tolerance ($\pm$ ml): 10.0	1 pc	2,062

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-10

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹550
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass			
GLT-GW03-049	Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.2 · Tolerance ($\pm$ ml): 0.2 · Height (mm): 160	1 pc (case of 2)	550
GLT-GW03-050	Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 25 · Graduation (ml): 0.5 · Tolerance ($\pm$ ml): 0.5 · Height (mm): 186	1 pc (case of 2)	619
GLT-GW03-051	Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 50 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 1.0 · Height (mm): 206	1 pc (case of 2)	688
GLT-GW03-052	Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 100 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 1.0 · Height (mm): 285	1 pc (case of 2)	819
GLT-GW03-053	Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 250 · Graduation (ml): 2.0 · Tolerance ($\pm$ ml): 2.0 · Height (mm): 370	1 pc	1,500

Full specifications and every pack size: Product Catalogue 2026-27, pages **485-491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW03-11

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass & Cylinders, graduated with round base and spout, amber stain graduation, Class 'A', with I/C stopper, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

4
PRODUCTS

4
PACK
OPTIONS

₹1,050
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass			
GLT-GW03-054	Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 500 · Graduation (ml): 5.0 · Tolerance ($\pm$ ml): 5.0 · Height (mm): 430	1 pc	1,994
GLT-GW03-055	Cylinders, graduated, with round base and spout, amber stain graduation, with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 1000 · Graduation (ml): 10.0 · Tolerance ($\pm$ ml): 10.0 · Height (mm): 500	1 pc	2,819
Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass			
GLT-GW03-056	Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.2 · Tolerance ($\pm$ ml): 0.1 · Height (mm): 160	1 pc (case of 2)	1,050
GLT-GW03-057	Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 25 · Graduation (ml): 0.5 · Tolerance ($\pm$ ml): 0.25 · Height (mm): 186	1 pc (case of 2)	1,188

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-12

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, graduated with round base and spout, amber stain graduation, Class 'A', with I/C stopper, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

4
PRODUCTS

4
PACK
OPTIONS

₹1,375
PRICE FROM

TYPICAL RESEARCH USE

Titrations

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass			
GLT-GW03-058	Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 50 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 0.5 · Height (mm): 206	1 pc (case of 2)	1,375
GLT-GW03-059	Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 100 · Graduation (ml): 1.0 · Tolerance ($\pm$ ml): 0.5 · Height (mm): 285	1 pc (case of 2)	1,638
GLT-GW03-060	Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 250 · Graduation (ml): 2.0 · Tolerance ($\pm$ ml): 1.0 · Height (mm): 370	1 pc	2,544
GLT-GW03-061	Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 500 · Graduation (ml): 5.0 · Tolerance ($\pm$ ml): 2.5 · Height (mm): 430	1 pc	3,300

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-13

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Cylinders, graduated with round base and spout, amber stain graduation, Class 'A', with I/C stopper, 3.3 borosilicate glass, Cylinders Crow with I/C Stopper, base and spout & more

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

6
PRODUCTS

6
PACK
OPTIONS

₹194
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass			
GLT-GW03-062	Cylinders, graduated with round base and spout, amber stain graduation, Class 'A' (with certificate), with I/C stopper, 3.3 borosilicate glass — Capacity (ml): 1000 · Graduation (ml): 10.0 · Tolerance ($\pm$ ml): 5.0 · Height (mm): 500	1 pc	3,919
Cylinders Crow with I/C Stopper, base and spout			
GLT-GW03-063	Cylinders Crow with I/C Stopper, base and spout — Capacity (ml): 50	1 pc (case of 2)	675
GLT-GW03-064	Cylinders Crow with I/C Stopper, base and spout — Capacity (ml): 100	1 pc (case of 2)	812
Pipettes, Milk testing, Bacteriological, graduated, 3.3 borosilicate glass			
GLT-GW03-065	Pipettes, Milk testing, Bacteriological, graduated, 3.3 borosilicate glass — Capacity (ml): 1.1 · Tolerance ($\pm$ ml): 0.025 · Marked at (ml): 0.5, 1.1	1 pc (case of 10)	194
GLT-GW03-066	Pipettes, Milk testing, Bacteriological, graduated, 3.3 borosilicate glass — Capacity (ml): 2.2 · Tolerance ($\pm$ ml): 0.040 · Marked at (ml): 1.0, 2.0,	1 pc (case of 10)	225
Pipettes, Milk testing Gerber, 3.3 borosilicate glass			
GLT-GW03-067	Pipettes, Milk testing Gerber, 3.3 borosilicate glass — Capacity (ml): 10.75 · Size (mm): 140	1 pc (case of 80)	194

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-14

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation, Class 'A

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹469
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation, Class 'A (with certificate)			
GLT-GW03-068	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation, Class 'A (with certificate) — Capacity (ml): 1 · Graduation (ml): 0.01 · Tolerance ($\pm$ ml): 0.006 · Colour code: Yellow	1 pc (case of 20)	469
GLT-GW03-069	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation, Class 'A (with certificate) — Capacity (ml): 2 · Graduation (ml): 0.02 · Tolerance ($\pm$ ml): 0.010 · Colour code: Black	1 pc (case of 20)	500
GLT-GW03-070	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation, Class 'A (with certificate) — Capacity (ml): 5 · Graduation (ml): 0.05 · Tolerance ($\pm$ ml): 0.030 · Colour code: Red	1 pc (case of 20)	575
GLT-GW03-071	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation, Class 'A (with certificate) — Capacity (ml): 10 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.050 · Colour code: Orange	1 pc (case of 20)	606
GLT-GW03-072	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation, Class 'A (with certificate) — Capacity (ml): 25 · Graduation (ml): 0.20 · Tolerance ($\pm$ ml): 0.100 · Colour code: White	1 pc (case of 20)	725

Full specifications and every pack size: Product Catalogue 2026-27, pages **485-491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW03-15

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹156
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation			
GLT-GW03-073	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation — Capacity (ml): 1 · Graduation (ml): 0.01 · Tolerance ($\pm$ ml): 0.01 · Colour code: Yellow	1 pc (case of 20)	156
GLT-GW03-074	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation — Capacity (ml): 2 · Graduation (ml): 0.02 · Tolerance ($\pm$ ml): 0.02 · Colour code: Black	1 pc (case of 20)	175
GLT-GW03-075	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation — Capacity (ml): 5 · Graduation (ml): 0.05 · Tolerance ($\pm$ ml): 0.05 · Colour code: Red	1 pc (case of 20)	188
GLT-GW03-076	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation — Capacity (ml): 10 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.10 · Colour code: Orange	1 pc (case of 20)	206
GLT-GW03-077	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation — Capacity (ml): 25 · Graduation (ml): 0.20 · Tolerance ($\pm$ ml): 0.20 · Colour code: White	1 pc (case of 20)	325

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-16

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Pipettes, Measuring, Mohr type, amber colour graduation, Class 'A & Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

6
PRODUCTS

6
PACK
OPTIONS

₹181
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Pipettes, Measuring, Mohr type, amber colour graduation, Class 'A (with certificate)			
GLT-GW03-078	Pipettes, Measuring, Mohr type, amber colour graduation, Class 'A (with certificate) — Capacity (ml): 1 · Subdivision (ml): 0.01 · Tolerance ($\pm$ ml): 0.006	1 pc (case of 20)	450
GLT-GW03-079	Pipettes, Measuring, Mohr type, amber colour graduation, Class 'A (with certificate) — Capacity (ml): 2 · Subdivision (ml): 0.02 · Tolerance ($\pm$ ml): 0.010	1 pc (case of 20)	475
GLT-GW03-080	Pipettes, Measuring, Mohr type, amber colour graduation, Class 'A (with certificate) — Capacity (ml): 5 · Subdivision (ml): 0.05 · Tolerance ($\pm$ ml): 0.030	1 pc (case of 20)	562
GLT-GW03-081	Pipettes, Measuring, Mohr type, amber colour graduation, Class 'A (with certificate) — Capacity (ml): 10 · Subdivision (ml): 0.10 · Tolerance ($\pm$ ml): 0.050	1 pc (case of 20)	600
Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass			
GLT-GW03-082	Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass — Capacity (ml): 0.1 · Graduation (ml): 0.01 · Tolerance ($\pm$ ml): 0.01	1 pc (case of 20)	181
GLT-GW03-083	Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass — Capacity (ml): 0.2 · Graduation (ml): 0.01 · Tolerance ($\pm$ ml): 0.01	1 pc (case of 20)	181

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-17

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass & Pipettes, Graduated, Serological, White colour graduation

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

7
PRODUCTS

7
PACK
OPTIONS

₹162
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass			
GLT-GW03-084	Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass — Capacity (ml): 1 · Graduation (ml): 0.01 · Tolerance ($\pm$ ml): 0.01	1 pc (case of 20)	162
GLT-GW03-085	Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass — Capacity (ml): 2 · Graduation (ml): 0.02 · Tolerance ($\pm$ ml): 0.02	1 pc (case of 20)	181
GLT-GW03-086	Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass — Capacity (ml): 5 · Graduation (ml): 0.05 · Tolerance ($\pm$ ml): 0.05	1 pc (case of 20)	188
GLT-GW03-087	Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 20)	200
GLT-GW03-088	Pipettes, Measuring, Mohr type, amber colour graduation, 3.3 borosilicate glass — Capacity (ml): 25 · Graduation (ml): 0.20 · Tolerance ($\pm$ ml): 0.20	1 pc (case of 12)	362
Pipettes, Graduated, Serological, White colour graduation			
GLT-GW03-089	Pipettes, Graduated, Serological, White colour graduation — Capacity (ml): 1 · Graduation (ml): 0.01 · Tolerance ($\pm$ ml): 0.01	1 pc (case of 20)	162
GLT-GW03-090	Pipettes, Graduated, Serological, White colour graduation — Capacity (ml): 2 · Graduation (ml): 0.02 · Tolerance ($\pm$ ml): 0.02	1 pc (case of 20)	175

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-18

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Pipettes, Graduated, Serological, White colour graduation & Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

5
PRODUCTS

5
PACK
OPTIONS

₹162
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Pipettes, Graduated, Serological, White colour graduation			
GLT-GW03-091	Pipettes, Graduated, Serological, White colour graduation — Capacity (ml): 5 · Graduation (ml): 0.05 · Tolerance ($\pm$ ml): 0.05	1 pc (case of 20)	194
GLT-GW03-092	Pipettes, Graduated, Serological, White colour graduation — Capacity (ml): 10 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.10	1 pc (case of 20)	206
Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation (upto tip), 3.3 borosilicate glass			
GLT-GW03-093	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation (upto tip), 3.3 borosilicate glass — Capacity (ml): 1 · Graduation (ml): 0.01 · Tolerance ($\pm$ ml): 0.01 · Colour code: Yellow	1 pc (case of 20)	162
GLT-GW03-094	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation (upto tip), 3.3 borosilicate glass — Capacity (ml): 2 · Graduation (ml): 0.02 · Tolerance ($\pm$ ml): 0.02 · Colour code: Black	1 pc (case of 20)	175
GLT-GW03-095	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation (upto tip), 3.3 borosilicate glass — Capacity (ml): 5 · Graduation (ml): 0.05 · Tolerance ($\pm$ ml): 0.05 · Colour code: Red	1 pc (case of 20)	194

Full specifications and every pack size: Product Catalogue 2026-27, pages **485-491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW03-19

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation, 3.3 borosilicate glass & Pipettes, Transfer, Volumetric, colour Coded, Class A

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

6
PRODUCTS

6
PACK
OPTIONS

₹200
PRICE FROM

TYPICAL RESEARCH USE

Titrations

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation (upto tip), 3.3 borosilicate glass			
GLT-GW03-096	Pipettes, Graduated, colour Coded, complete delivery, amber colour graduation (upto tip), 3.3 borosilicate glass — Capacity (ml): 10 · Graduation (ml): 0.10 · Tolerance ($\pm$ ml): 0.10 · Colour code: Orange	1 pc (case of 20)	200
Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate)			
GLT-GW03-097	Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate) — Capacity (ml): 1 · Tolerance ($\pm$ ml): 0.007 · Colour code: Blue	1 pc (case of 20)	350
GLT-GW03-098	Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate) — Capacity (ml): 2 · Tolerance ($\pm$ ml): 0.010 · Colour code: Orange	1 pc (case of 20)	350
GLT-GW03-099	Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate) — Capacity (ml): 5 · Tolerance ($\pm$ ml): 0.015 · Colour code: White	1 pc (case of 20)	400
GLT-GW03-100	Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate) — Capacity (ml): 10 · Tolerance ($\pm$ ml): 0.020 · Colour code: Red	1 pc (case of 20)	475
GLT-GW03-101	Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate) — Capacity (ml): 20 · Tolerance ($\pm$ ml): 0.030 · Colour code: Yellow	1 pc (case of 20)	606

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-20

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Pipettes, Transfer, Volumetric, colour Coded, Class A & Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

7
PRODUCTS

7
PACK
OPTIONS

₹162
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate)			
GLT-GW03-102	Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate) — Capacity (ml): 25 · Tolerance ($\pm$ ml): 0.030 · Colour code: Blue	1 pc (case of 20)	662
GLT-GW03-103	Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate) — Capacity (ml): 50 · Tolerance ($\pm$ ml): 0.050 · Colour code: Red	1 pc (case of 20)	744
GLT-GW03-104	Pipettes, Transfer, Volumetric, colour Coded, Class A (with y certificate) — Capacity (ml): 100 · Tolerance ($\pm$ ml): 0.080 · Colour code: Yellow	1 pc (case of 20)	1,131
Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass			
GLT-GW03-105	Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass — Capacity (ml): 1 · Tolerance ($\pm$ ml): 0.015 · Colour code: Blue	1 pc (case of 20)	162
GLT-GW03-106	Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass — Capacity (ml): 2 · Tolerance ($\pm$ ml): 0.020 · Colour code: Orange	1 pc (case of 20)	162
GLT-GW03-107	Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass — Capacity (ml): 5 · Tolerance ($\pm$ ml): 0.030 · Colour code: White	1 pc (case of 20)	169
GLT-GW03-108	Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass — Capacity (ml): 10 · Tolerance ($\pm$ ml): 0.040 · Colour code: Red	1 pc (case of 20)	206

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-21

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW03 BURETTES, PIPETTES & MEASURING CYLINDERS

Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass & Pipettes, Pasteur, 3.3 borosilicate glass

Illustrative image

Burettes, pipettes and measuring cylinders, including Class A items listed with certificate.

6
PRODUCTS

6
PACK
OPTIONS

₹31
PRICE FROM

TYPICAL RESEARCH USE

Titration

Accurate liquid transfer and dilution

Measuring liquid volumes

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass			
GLT-GW03-109	Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass — Capacity (ml): 20 · Tolerance ($\pm$ ml): 0.060 · Colour code: Yellow	1 pc (case of 12)	238
GLT-GW03-110	Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass — Capacity (ml): 25 · Tolerance ($\pm$ ml): 0.060 · Colour code: Blue	1 pc (case of 12)	250
GLT-GW03-111	Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass — Capacity (ml): 50 · Tolerance ($\pm$ ml): 0.100 · Colour code: Red	1 pc (case of 12)	306
GLT-GW03-112	Pipettes, Transfer, Volumetric, colour Coded, 3.3 borosilicate glass — Capacity (ml): 100 · Tolerance ($\pm$ ml): 0.160 · Colour code: Yellow	1 pc (case of 12)	444
Pipettes, Pasteur, 3.3 borosilicate glass			
GLT-GW03-113	Pipettes, Pasteur, 3.3 borosilicate glass — Size ; (mm): 140	1 pc (case of 100)	31
GLT-GW03-114	Pipettes, Pasteur, 3.3 borosilicate glass — Size ; (mm): 230	1 pc (case of 100)	38

Full specifications and every pack size: Product Catalogue 2026–27, pages **485–491**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW03-22

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113