

GW08 GLASS VIALS, SLIDES, CUVETTES & CONSUMABLES

Glass Vial / Sample Vial, Glass Vial / Culture Tube & more

Illustrative image

Glass vials, slides and cover glasses, optical cells, capillaries, ampoules and counting chambers.

13
PRODUCTS

16
PACK
OPTIONS

₹625
PRICE FROM

TYPICAL RESEARCH USE

Microscopy

Spectrophotometry

Sample storage and cell counting

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Glass Vial / Sample Vial			
GLT-GW08-001 +1	Glass Vial / Sample Vial — Capacity: 2ml	Pack of 500 vials · Pack of 100 vials	from 625
GLT-GW08-002 +1	Glass Vial / Sample Vial — Capacity: 5ml	Pack of 500 vials · Pack of 100 vials	from 1,125
GLT-GW08-003 +1	Glass Vial / Sample Vial — Capacity: 10ml	Pack of 500 vials · Pack of 100 vials	from 1,375
GLT-GW08-007	Glass Vial / Sample Vial — Capacity: 50ml	Pack of 100 vials	3,750
Glass Vial / Sample Vial (Amber)			
GLT-GW08-008	Glass Vial / Sample Vial (Amber) — Capacity: 5ml	Pack of 100 vials	1,250
GLT-GW08-009	Glass Vial / Sample Vial (Amber) — Capacity: 10ml	Pack of 100 vials	1,500
Glass Vial / Culture Tube			
GLT-GW08-010	Glass Vial / Culture Tube — Capacity: 5ml	Pack of 100 vials	2,750
GLT-GW08-011	Glass Vial / Culture Tube — Capacity: 15ml	Pack of 100 vials	4,000
GLT-GW08-012	Glass Vial / Culture Tube — Capacity: 30ml	Pack of 100 vials	4,250
GLT-GW08-013	Glass Vial / Culture Tube — Capacity: 50ml	Pack of 100 vials	4,750
Capillary			
GLT-GW08-014	Capillary — Type: (pack of 100 capillaries)	Pack (see description)	1,000
GLT-GW08-015	Capillary — Type: (Vial of 100 Capillary)	Pack (see description)	1,750
GLT-GW08-016	Capillary — Type: Melting Point	Pack (see description)	2,625

Full specifications and every pack size: Product Catalogue 2026–27, pages **509–513**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW08-01

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW08 GLASS VIALS, SLIDES, CUVETTES & CONSUMABLES

Capillary & Glass Syringe

Illustrative image

Glass vials, slides and cover glasses, optical cells, capillaries, ampoules and counting chambers.

13
PRODUCTS

13
PACK
OPTIONS

₹1,375
PRICE FROM

TYPICAL RESEARCH USE

Microscopy

Spectrophotometry

Sample storage and cell counting

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Capillary			
GLT-GW08-017	Capillary — Type: Vial/100 Capillary) · option 1	Pack (see description)	2,625
GLT-GW08-018	Capillary — Type: Vial/100 Capillary) · option 2	Pack (see description)	2,625
GLT-GW08-019	Capillary — Type: Vial/100 Capillary)	Pack (see description)	2,625
Glass Syringe			
GLT-GW08-020	Glass Syringe — (Luer Mount) "INSULIN" · Capacity: 1 ml	Pack of 10	1,375
GLT-GW08-021	Glass Syringe — (Luer Mount) "TUBERCULIN" · Capacity: 1 ml	Pack of 10	1,750
GLT-GW08-022	Glass Syringe — Luer Mount · Capacity: 2 ml	Pack of 10	1,500
GLT-GW08-023	Glass Syringe — Luer Mount · Capacity: 5 ml	Pack of 10	1,938
GLT-GW08-024	Glass Syringe — Luer Mount · Capacity: 10 ml	Pack of 10	2,250
GLT-GW08-025	Glass Syringe — Luer Mount · Capacity: 20 ml	Pack of 6	1,800
GLT-GW08-026	Glass Syringe — Luer Lock · Capacity: 2 ml	Pack of 10	1,500
GLT-GW08-027	Glass Syringe — Luer Lock · Capacity: 5 ml	Pack of 10	1,875
GLT-GW08-028	Glass Syringe — Luer Lock · Capacity: 10 ml	Pack of 10	2,625
GLT-GW08-029	Glass Syringe — Luer Lock · Capacity: 20 ml	Pack of 6	2,700

Full specifications and every pack size: Product Catalogue 2026–27, pages **509–513**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW08-02

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW08 GLASS VIALS, SLIDES, CUVETTES & CONSUMABLES

Optical Cells

Illustrative image

Glass vials, slides and cover glasses, optical cells, capillaries, ampoules and counting chambers.

8
PRODUCTS

8
PACK
OPTIONS

₹8,000
PRICE FROM

TYPICAL RESEARCH USE

Microscopy

Spectrophotometry

Sample storage and cell counting

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Optical Cells (Cuvettes)			
GLT-GW08-030	Optical Cells (Cuvettes) — Optical path Length (mm): 1mm · External Size (mm): 12.5×3.5×45 · Capacity (ml): 0.35ml	1 matched pair	25,781
GLT-GW08-031	Optical Cells (Cuvettes) — Optical path Length (mm): 10mm · External Size (mm): 12.5×12.5×45 · Capacity (ml): 0.7ml	1 matched pair	21,000
GLT-GW08-032	Optical Cells (Cuvettes) — Optical path Length (mm): (Black)10mm Oblique midget · External Size (mm): 12.5×12.5×45 · Capacity (ml): 1.00ml	1 matched pair	19,500
GLT-GW08-033	Optical Cells (Cuvettes) — Optical path Length (mm): Oblique midget · External Size (mm): 12.5×12.5×45 · Capacity (ml): 1.40ml	1 matched pair	19,000
GLT-GW08-034	Optical Cells (Cuvettes) — Optical path Length (mm): 10mm · External Size (mm): 12.5×12.5×45 · Capacity (ml): 3.50ml	1 matched pair	8,000
GLT-GW08-035	Optical Cells (Cuvettes) — Optical path Length (mm): 20mm · External Size (mm): 12.5×22.5×45 · Capacity (ml): 7.00ml	1 matched pair	16,000
GLT-GW08-036	Optical Cells (Cuvettes) — Optical path Length (mm): 40mm · External Size (mm): 12.5×42.5×45 · Capacity (ml): 14.00ml	1 matched pair	19,000
GLT-GW08-037	Optical Cells (Cuvettes) — Optical path Length (mm): 50mm · External Size (mm): 12.5×52.5×45 · Capacity (ml): 17.50ml	1 matched pair	22,000

Full specifications and every pack size: Product Catalogue 2026–27, pages **509–513**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW08-03

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW08 GLASS VIALS, SLIDES, CUVETTES & CONSUMABLES

Optical Cells

Illustrative image

Glass vials, slides and cover glasses, optical cells, capillaries, ampoules and counting chambers.

7
PRODUCTS

7
PACK
OPTIONS

₹4,500
PRICE FROM

TYPICAL RESEARCH USE

Microscopy

Spectrophotometry

Sample storage and cell counting

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Optical Cells (Cuvettes)			
GLT-GW08-038	Optical Cells (Cuvettes) — Optical path Length (mm): 100mm · External Size (mm): 12.5×102.5×45 · Capacity (ml): 35.00ml	1 matched pair	41,000
Optical Cells			
GLT-GW08-039	Optical Cells — Optical path Length (mm): Oblique midget10mm · External Size (mm): 12.5×12.5×45 · Capacity (ml): 0.7ml	1 matched pair	16,312
GLT-GW08-040	Optical Cells — Optical path Length (mm): " " (Black)10mm · External Size (mm): 12.5×12.5×45 · Capacity (ml): 1.40ml	1 matched pair	15,625
GLT-GW08-041	Optical Cells — Optical path Length (mm): 10mm · External Size (mm): 12.5×12.5×45 · Capacity (ml): 3.50ml	1 matched pair	4,500
GLT-GW08-042	Optical Cells — Optical path Length (mm): 20mm · External Size (mm): 12.5×22.5×45 · Capacity (ml): 7.00ml	1 matched pair	6,750
GLT-GW08-043	Optical Cells — Optical path Length (mm): 40mm · External Size (mm): 12.5×42.5×45 · Capacity (ml): 14.00ml	1 matched pair	15,500
GLT-GW08-044	Optical Cells — Optical path Length (mm): 50mm · External Size (mm): 12.5×52.5×45 · Capacity (ml): 17.50ml	1 matched pair	17,625

Full specifications and every pack size: Product Catalogue 2026–27, pages **509–513**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW08-04

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW08 GLASS VIALS, SLIDES, CUVETTES & CONSUMABLES

Optical Cells, Cells For Fluorimeter-Quartz & more

Illustrative image

Glass vials, slides and cover glasses, optical cells, capillaries, ampoules and counting chambers.

9
PRODUCTS

9
PACK
OPTIONS

₹90
PRICE FROM

TYPICAL RESEARCH USE

Microscopy

Spectrophotometry

Sample storage and cell counting

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Optical Cells			
GLT-GW08-045	Optical Cells — Optical path Length (mm): 100mm · External Size (mm): 12.5×102.5×45 · Capacity (ml): 35.00ml	1 matched pair	21,500
Cells For Fluorimeter- Quartz			
GLT-GW08-046	Cells For Fluorimeter- Quartz — path Length (mm): 10mm · Type: Macro · Capacity (ml): 3.50ml	1 matched pair	26,312
GLT-GW08-047	Cells For Fluorimeter- Quartz — path Length (mm): 10mm · Type: Macro · Capacity (ml): 104ml	1 matched pair	40,625
Neubauer Counting Chamber			
GLT-GW08-048	Neubauer Counting Chamber — option 1	Pack of 2	4,050
GLT-GW08-049	Neubauer Counting Chamber — option 2	Pack of 2	4,050
Lyophil Vial/Ampules Glass			
GLT-GW08-050	Lyophil Vial/Ampules Glass — Capacity: 2ml · 10×122×7mm	Pack of 100 ampul	3,000
GLT-GW08-051	Lyophil Vial/Ampules Glass — Capacity: 5ml · 14×122×7mm	Pack of 100 ampul	3,000
Microscopic Glass Slides			
GLT-GW08-052	Microscopic Glass Slides — Plane Glass Slide (PIC-2) · Size; Lx W x thickness: 75×25×1.45mm	1 pc	90
GLT-GW08-053	Microscopic Glass Slides — Plane Glass Slide (PIC-1) · Size; Lx W x thickness: 75×25×1.35mm	1 pc	100

Full specifications and every pack size: Product Catalogue 2026-27, pages 509-513.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-GW08-05

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW08 GLASS VIALS, SLIDES, CUVETTES & CONSUMABLES

Microscopic Glass Slides & Microscopic Cover Glass

Illustrative image

Glass vials, slides and cover glasses, optical cells, capillaries, ampoules and counting chambers.

10
PRODUCTS

10
PACK
OPTIONS

₹120
PRICE FROM

TYPICAL RESEARCH USE

Microscopy

Spectrophotometry

Sample storage and cell counting

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Microscopic Glass Slides			
GLT-GW08-054	Microscopic Glass Slides — Plane Glass Slide (Deluxe) · Size; Lx W x thickness: 75×25×1.25mm	1 pc	120
GLT-GW08-055	Microscopic Glass Slides — Plane Glass Slide (S.Deluxe) · Size; Lx W x thickness: 75×25×1.10mm	1 pc	150
GLT-GW08-056	Microscopic Glass Slides — Plane Glass Slide (S.S.Deluxe) · Size; Lx W x thickness: 75×25×1mm	1 pc	194
GLT-GW08-057	Microscopic Glass Slides — Forested Micro Slide · Size; Lx W x thickness: 75×25×1.35mm	1 pc	120
GLT-GW08-058	Microscopic Glass Slides — Forested Micro Slide · Size; Lx W x thickness: 75×25×1.25mm	1 pc	138
GLT-GW08-059	Microscopic Glass Slides — Forested Micro Slide · Size; Lx W x thickness: 75×25×1.10mm	1 pc	165
Microscopic Cover Glass(Cover Slips)			
GLT-GW08-060	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Square:18×18mm	10 g pack	200
GLT-GW08-061	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Square: 22×22mm	10 g pack	200
GLT-GW08-062	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Rectangular:22×30mm	10 g pack	225
GLT-GW08-063	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Rectangular:22×40mm	10 g pack	225

Full specifications and every pack size: Product Catalogue 2026–27, pages **509–513**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW08-06

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW08 GLASS VIALS, SLIDES, CUVETTES & CONSUMABLES

Microscopic Cover Glass & Cap,

Illustrative image

Glass vials, slides and cover glasses, optical cells, capillaries, ampoules and counting chambers.

10
PRODUCTS

10
PACK
OPTIONS

₹225
PRICE FROM

TYPICAL RESEARCH USE

Microscopy

Spectrophotometry

Sample storage and cell counting

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Microscopic Cover Glass(Cover Slips)			
GLT-GW08-064	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Rectangular:22×50mm	10 g pack	225
GLT-GW08-065	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Rectangular:22×60mm	10 g pack	262
GLT-GW08-066	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Rectangular:24×50mm	10 g pack	262
GLT-GW08-067	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Rectangular:24×60mm	10 g pack	262
GLT-GW08-068	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Circular: 18mm	10 g pack	350
GLT-GW08-069	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Circular: 22mm	10 g pack	350
GLT-GW08-070	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Circular: 16mm	10 g pack	362
GLT-GW08-071	Microscopic Cover Glass(Cover Slips) — Micro Cover Glass · Sizes: Circular: 12mm	10 g pack	362
GLT-GW08-072	Microscopic Cover Glass(Cover Slips) — Haemocytometer Cover Glass · Sizes: 20mmx25mm	10 g pack	300
Cap, (Spare Cap for Glass Vial)			
GLT-GW08-073	Cap, (Spare Cap for Glass Vial) — Capacity For Vial: 5ml · Fitting tube dia: 16 mm	Pack of 100	1,250

Full specifications and every pack size: Product Catalogue 2026–27, pages **509–513**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW08-07

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW08 GLASS VIALS, SLIDES, CUVETTES & CONSUMABLES

Cap, & Bottles, relative density, with capillary bore, PTFE stopper, 3.3 borosilicate glass

Illustrative image

Glass vials, slides and cover glasses, optical cells, capillaries, ampoules and counting chambers.

8
PRODUCTS

8
PACK
OPTIONS

₹581
PRICE FROM

TYPICAL RESEARCH USE

Microscopy

Spectrophotometry

Sample storage and cell counting

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cap, (Spare Cap for Glass Vial)			
GLT-GW08-074	Cap, (Spare Cap for Glass Vial) — Capacity For Vial: 15ml · Fitting tube dia: 25mm	Pack of 100	1,500
GLT-GW08-075	Cap, (Spare Cap for Glass Vial) — Capacity For Vial: 30ml · Fitting tube dia: 25mm	Pack of 100	1,500
GLT-GW08-076	Cap, (Spare Cap for Glass Vial) — Capacity For Vial: 50ml · Fitting tube dia: 25mm	Pack of 100	1,500
GLT-GW08-077	Cap, (Spare Cap for Glass Vial) — Capacity For Vial: 5ml · Fitting tube dia: 20mm	Pack of 1000	1,500
GLT-GW08-078	Cap, (Spare Cap for Glass Vial) — Capacity For Vial: 2ml · Fitting tube dia: 13mm	Pack of 1000	1,250
GLT-GW08-079	Cap, (Spare Cap for Glass Vial) — Capacity For Vial: 5ml & 10ml · Fitting tube dia: 20mm	Pack of 1000	2,250
Bottles, relative density, with capillary bore, PTFE stopper, 3.3 borosilicate glass			
GLT-GW08-080	Bottles, relative density, with capillary bore, PTFE stopper, 3.3 borosilicate glass — 10 · Tolerance ($\pm$ ml): 1 · Dimensions (mm): 24 × 46 · Stopper size: 10/15	1 pc (case of 2)	581
GLT-GW08-081	Bottles, relative density, with capillary bore, PTFE stopper, 3.3 borosilicate glass — 25 · Tolerance ($\pm$ ml): 2 · Dimensions (mm): 38 × 56 · Stopper size: 10/15	1 pc (case of 2)	631

Full specifications and every pack size: Product Catalogue 2026–27, pages **509–513**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW08-08

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

GW08 GLASS VIALS, SLIDES, CUVETTES & CONSUMABLES

Bottles, relative density, with capillary bore, PTFE stopper, 3.3 borosilicate glass

Illustrative image

Glass vials, slides and cover glasses, optical cells, capillaries, ampoules and counting chambers.

1
PRODUCTS

1
PACK
OPTIONS

₹688
PRICE FROM

TYPICAL RESEARCH USE

Microscopy

Spectrophotometry

Sample storage and cell counting

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Bottles, relative density, with capillary bore, PTFE stopper, 3.3 borosilicate glass			
GLT-GW08-082	Bottles, relative density, with capillary bore, PTFE stopper, 3.3 borosilicate glass — 50 · Tolerance ($\pm$ ml): 3 · Dimensions (mm): 49 × 66 · Stopper size: 10/15	1 pc (case of 2)	688

Full specifications and every pack size: Product Catalogue 2026–27, pages **509–513**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-GW08-09

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113