

IN03 POWER SUPPLIES

Analog Model: Fixed Voltage, Digital PS 100 & more

Illustrative image

Power supplies for electrophoresis and electrotransfer with constant-voltage and constant-current modes.

8
PRODUCTS

8
PACK
OPTIONS

₹7,120
PRICE FROM

TYPICAL RESEARCH USE

Powering agarose and PAGE runs

Powering blot transfers

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Analog Model: Fixed Voltage			
GLT-IN03-001	PS 100: Mini Model 50, 100, 150, 200, 250 V Fixed 100 mA with on/off LED Indicator	1 EA	7,120
GLT-IN03-002	PS 200: Mini Model 100, 150, 200, 250 V Fixed, 200 mA	1 EA	22,740
GLT-IN03-003	Mini Powerpack with timer	1 EA	27,960
Digital PS 100			
GLT-IN03-004	Digital PS 100: EPS 300 Model ; 0-300 V; 0-100 mA ; Variable with timer Constant Current ; Constant Voltage ; Digital display	1 EA	30,410
Digital PS 500			
GLT-IN03-005	EPS 500 Model ;0-500 V; 0-500 mA; Variable with timer Constant Current ; Constant Voltage ; Digital display	1 EA	60,970
GLT-IN03-006	Digital PS 3000: Electrotransfer Model ; 10-250 V ; 3000 mA; Variable with timer Constant Voltage ; Digital display	1 EA	73,070
GLT-IN03-007	Euro Model Digital Power supply light weight; 10-300 V; 4-400 mA Max 75W, 4 outlets; Variable with timer; Constant Current; Constant Voltage; Digital display	1 EA	72,350
GLT-IN03-008	300 V/150mA Power Supply	1 EA	47,450

Full specifications and every pack size: Product Catalogue 2026-27, pages 96.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-IN03-01

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113