

LW02 STAINLESS STEEL & METAL LABWARE

Enamel Trays, Trays S.S & more

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

13
PRODUCTS

13
PACK
OPTIONS

₹550
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Enamel Trays			
GLT-LW02-001	Enamel Trays — Size (WXL): 6×8"	Pack of 2	1,050
GLT-LW02-002	Enamel Trays — Size (WXL): 8×10"	Pack of 2	1,350
GLT-LW02-003	Enamel Trays — Size (WXL): 10×12"	Pack of 2	1,700
GLT-LW02-004	Enamel Trays — Size (WXL): 12×15"	Pack of 2	2,300
GLT-LW02-005	Enamel Trays — Size (WXL): 12×18"	Pack of 2	3,100
GLT-LW02-006	Enamel Trays — Size (WXL): 18×24"	Pack of 2	4,100
Trays S.S			
GLT-LW02-007	Trays S.S — Size (WXL): 6×8"	Pack of 2	1,100
GLT-LW02-008	Trays S.S — Size (WXL): 8×10"	Pack of 2	1,550
GLT-LW02-009	Trays S.S — Size (WXL): 10×12"	Pack of 2	1,900
GLT-LW02-010	Trays S.S — Size (WXL): 10×14"	Pack of 2	2,300
GLT-LW02-011	Trays S.S — Size (WXL): 12×15"	Pack of 2	3,200
S.S. Test Tube Stand			
GLT-LW02-012	S.S. Test Tube Stand — Size (Dia x Hole): 12 × 12	Pack of 2	550
GLT-LW02-013	S.S. Test Tube Stand — Size (Dia x Hole): 12 × 24	Pack of 2	812

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-01

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

S.S. Test Tube Stand

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

14
PRODUCTS

14
PACK
OPTIONS

₹650
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
S.S. Test Tube Stand			
GLT-LW02-014	S.S. Test Tube Stand — Size (Dia x Hole): 12 × 36	Pack of 2	975
GLT-LW02-015	S.S. Test Tube Stand — Size (Dia x Hole): 12 × 48	Pack of 2	1,125
GLT-LW02-016	S.S. Test Tube Stand — Size (Dia x Hole): 15 × 12	Pack of 2	650
GLT-LW02-017	S.S. Test Tube Stand — Size (Dia x Hole): 15 × 24	Pack of 2	875
GLT-LW02-018	S.S. Test Tube Stand — Size (Dia x Hole): 15 × 36	Pack of 2	1,000
GLT-LW02-019	S.S. Test Tube Stand — Size (Dia x Hole): 15 × 48	Pack of 2	1,200
GLT-LW02-020	S.S. Test Tube Stand — Size (Dia x Hole): 18 × 12	Pack of 2	750
GLT-LW02-021	S.S. Test Tube Stand — Size (Dia x Hole): 18 × 24	Pack of 2	1,062
GLT-LW02-022	S.S. Test Tube Stand — Size (Dia x Hole): 18 × 36	Pack of 2	1,200
GLT-LW02-023	S.S. Test Tube Stand — Size (Dia x Hole): 18 × 48	Pack of 2	1,375
GLT-LW02-024	S.S. Test Tube Stand — Size (Dia x Hole): 25 × 12	Pack of 2	875
GLT-LW02-025	S.S. Test Tube Stand — Size (Dia x Hole): 25 × 24	Pack of 2	1,200
GLT-LW02-026	S.S. Test Tube Stand — Size (Dia x Hole): 25 × 36	Pack of 2	1,600
GLT-LW02-027	S.S. Test Tube Stand — Size (Dia x Hole): 25 × 48	Pack of 2	2,062

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-02

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Kidney Trays S.S., Scoops S.S. & more

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

12
PRODUCTS

12
PACK
OPTIONS

₹938
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Kidney Trays S.S.			
GLT-LW02-028	Kidney Trays S.S. — Size (Length): 6"	Pack of 6	938
GLT-LW02-029	Kidney Trays S.S. — Size (Length): 8"	Pack of 6	1,088
GLT-LW02-030	Kidney Trays S.S. — Size (Length): 10"	Pack of 6	1,275
GLT-LW02-031	Kidney Trays S.S. — Size (Length): 12"	Pack of 6	1,575
Scoops S.S.			
GLT-LW02-032	Scoops S.S. — Capacity: 250gm	Pack of 2	2,550
GLT-LW02-033	Scoops S.S. — Capacity: ½ kg	Pack of 2	2,700
GLT-LW02-034	Scoops S.S. — Capacity: 1kg	Pack of 2	3,200
GLT-LW02-035	Scoops S.S. — Capacity: 1½ kg	Pack of 2	3,900
GLT-LW02-036	Scoops S.S. — Capacity: 2kg	Pack of 2	5,025
GLT-LW02-037	Scoops S.S. — Capacity: 2½ kg	Pack of 2	5,650
Vial Sealer / 'D' Sealer			
GLT-LW02-038	Vial Sealer / 'D' Sealer — Vial Neck Size: 13mm · Vial Sealer	1 pc	4,500
GLT-LW02-039	Vial Sealer / 'D' Sealer — Vial Neck Size: 20mm · Vial Sealer	1 pc	4,500

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-03

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Vial Sealer / 'D' Sealer & Test Tube Cap S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

10
PRODUCTS

10
PACK
OPTIONS

₹3,600
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Vial Sealer / 'D' Sealer			
GLT-LW02-040	Vial Sealer / 'D' Sealer — Vial Neck Size: 13mm · Vial DSealer	1 pc	4,500
GLT-LW02-041	Vial Sealer / 'D' Sealer — Vial Neck Size: 20mm · Vial DSealer	1 pc	4,500
GLT-LW02-042	Vial Sealer / 'D' Sealer — Vial Neck Size: 13mm · Vial Sealer Imp.	1 pc	26,250
GLT-LW02-043	Vial Sealer / 'D' Sealer — Vial Neck Size: 13mm · Vial DSealer Imp.	1 pc	26,250
GLT-LW02-044	Vial Sealer / 'D' Sealer — Vial Neck Size: 20mm · Vial Sealer Imp.	1 pc	26,250
GLT-LW02-045	Vial Sealer / 'D' Sealer — Vial Neck Size: 20mm · Vial DSealer Imp.	1 pc	26,250
Test Tube Cap S.S.			
GLT-LW02-046	Test Tube Cap S.S. — For test tube: 12×100mm	Pack of 24	3,600
GLT-LW02-047	Test Tube Cap S.S. — For test tube: 15×125mm	Pack of 24	4,500
GLT-LW02-048	Test Tube Cap S.S. — For test tube: 18×150mm	Pack of 24	5,700
GLT-LW02-049	Test Tube Cap S.S. — For test tube: 25×150mm	Pack of 24	6,600

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-04

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Pipette Container S.S., Laboratory Jack & more

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

9
PRODUCTS

9
PACK
OPTIONS

₹650
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Pipette Container S.S.			
GLT-LW02-050	Pipette Container S.S. — Size (Length x Dial: 16×2.5"	Pack of 2	2,150
GLT-LW02-051	Pipette Container S.S. — Size (Length x Dial: 18×3"	Pack of 2	2,300
Laboratory Jack			
GLT-LW02-052	Laboratory Jack — 100mm x 150mm	1 pc	4,525
GLT-LW02-053	Laboratory Jack — 150mm x 200mm	1 pc	5,575
GLT-LW02-054	Laboratory Jack — 200mm x 200mm	1 pc	6,600
Slide Carrier/Slide Trays			
GLT-LW02-055	Slide Carrier/Slide Trays — Description & No. of Place: Slide Carrier · Alu/S.S: Alum...	Pack of 2	650
GLT-LW02-056	Slide Carrier/Slide Trays — Description & No. of Place: " " · Alu/S.S: S.S. · option 1	Pack of 2	1,350
GLT-LW02-057	Slide Carrier/Slide Trays — Description & No. of Place: Slide Staining · Alu/S.S: Alum...	Pack of 2	950
GLT-LW02-058	Slide Carrier/Slide Trays — Description & No. of Place: " " · Alu/S.S: S.S. · option 2	Pack of 2	1,700

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-05

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Slide Staining Box & Forceps/Bull Dog Clamps/Bone Cutter

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

11
PRODUCTS

11
PACK
OPTIONS

₹1,375
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Slide Staining Box			
GLT-LW02-059	Slide Staining Box — Size (mm) (LxWx H): 130×90×40mm	Pack of 2	3,650
Forceps/Bull Dog Clamps/Bone Cutter			
GLT-LW02-060	Forceps/Bull Dog Clamps/Bone Cutter — Forceps 4"	Pack of 10	1,375
GLT-LW02-061	Forceps/Bull Dog Clamps/Bone Cutter — Forceps 5"	Pack of 10	1,500
GLT-LW02-062	Forceps/Bull Dog Clamps/Bone Cutter — Forceps 6"	Pack of 10	1,750
GLT-LW02-063	Forceps/Bull Dog Clamps/Bone Cutter — Forceps 8"	Pack of 10	2,000
GLT-LW02-064	Forceps/Bull Dog Clamps/Bone Cutter — Forceps 10"	Pack of 10	2,250
GLT-LW02-065	Forceps/Bull Dog Clamps/Bone Cutter — Forceps 12"	Pack of 10	2,500
GLT-LW02-066	Forceps/Bull Dog Clamps/Bone Cutter — Artery Forceps 6"	Pack of 10	3,250
GLT-LW02-067	Forceps/Bull Dog Clamps/Bone Cutter — Capsula Forceps 4"	Pack of 4	2,400
GLT-LW02-068	Forceps/Bull Dog Clamps/Bone Cutter — Bull dog Clamps (S)	Pack of 6	2,700
GLT-LW02-069	Forceps/Bull Dog Clamps/Bone Cutter — Bull dog clamps (M)	Pack of 6	2,700

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-06

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Forceps/Bull Dog Clamps/Bone Cutter, B.P Handle S.S. Deluxe Quality & more

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

10
PRODUCTS

10
PACK
OPTIONS

₹1,100
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Forceps/Bull Dog Clamps/Bone Cutter			
GLT-LW02-070	Forceps/Bull Dog Clamps/Bone Cutter — Bone Cutter 6"	Pack of 2	2,550
GLT-LW02-071	Forceps/Bull Dog Clamps/Bone Cutter — Bone cutter 8"	Pack of 2	3,150
B.P Handle S.S. Deluxe Quality			
GLT-LW02-072	B.P Handle S.S. Deluxe Quality — Number: 4 · Suitable for	Pack of 10	1,125
GLT-LW02-073	B.P Handle S.S. Deluxe Quality — Number: 3 · Suitable for	Pack of 10	1,125
Scissor S.S.			
GLT-LW02-074	Scissor S.S. — 4"	Pack of 10	1,750
GLT-LW02-075	Scissor S.S. — 5"	Pack of 10	2,000
GLT-LW02-076	Scissor S.S. — 6"	Pack of 10	2,250
GLT-LW02-077	Scissor S.S. — 8"	Pack of 4	1,100
GLT-LW02-078	Scissor S.S. — 3" Iris Scissor Curved	Pack of 4	3,100
GLT-LW02-079	Scissor S.S. — 3" Iris Scissor Straight	Pack of 4	3,100

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-07

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Spatula S.S & Tong S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

14
PRODUCTS

14
PACK
OPTIONS

₹700
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Spatula S.S			
GLT-LW02-080	Spatula S.S — 6" · Type: Micro Wire Type	Pack of 10	1,000
GLT-LW02-081	Spatula S.S — 6" · Type: Micro One	Pack of 10	1,000
GLT-LW02-082	Spatula S.S — 8" · Type: Wire Type	Pack of 10	1,062
GLT-LW02-083	Spatula S.S — 10" · Type: Wire Type	Pack of 10	1,375
GLT-LW02-084	Spatula S.S — 12" · Type: Wire Type	Pack of 10	1,562
GLT-LW02-085	Spatula S.S — 16" · Type: Wire Type	Pack of 10	2,000
GLT-LW02-086	Spatula S.S — 4" · Type: Spoon Type	Pack of 10	938
GLT-LW02-087	Spatula S.S — 6" · Type: Spoon Type	Pack of 10	1,125
GLT-LW02-088	Spatula S.S — 8" · Type: Spoon Type	Pack of 10	1,375
GLT-LW02-089	Spatula S.S — 10" · Type: Spoon Type	Pack of 10	1,56,250
GLT-LW02-090	Spatula S.S — 12" · Type: Spoon Type	Pack of 10	2,000
Tong S.S. (Heavy Duty)			
GLT-LW02-091	Tong S.S. (Heavy Duty) — 10'	Pack of 2	700
GLT-LW02-092	Tong S.S. (Heavy Duty) — 12'	Pack of 2	800
GLT-LW02-093	Tong S.S. (Heavy Duty) — 16'	Pack of 2	1,150

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-08

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Tong S.S., Stand & more

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

9
PRODUCTS

9
PACK
OPTIONS

₹1,550
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Tong S.S. (Heavy Duty)			
GLT-LW02-094	Tong S.S. (Heavy Duty) — 18"	Pack of 2	1,550
GLT-LW02-095	Tong S.S. (Heavy Duty) — 20"	Pack of 2	2,050
Stand (Retort)			
GLT-LW02-096	Stand (Retort) — 8×5×24"x12mm	Pack of 2	1,550
GLT-LW02-097	Stand (Retort) — 9×6×30"x12mm	Pack of 2	1,700
GLT-LW02-098	Stand (Retort) — 12×7×36"x12mm	Pack of 2	1,900
Laboratory Clamps			
GLT-LW02-099	Laboratory Clamps — Size (Aprox.): 14" · 3Finger with	Pack of 10 sets	5,250
GLT-LW02-100	Laboratory Clamps — Size (Aprox.): 14" · 4Finger with	Pack of 10 sets	5,250
GLT-LW02-101	Laboratory Clamps — Size (Aprox.): For 3/4"(19mm) · Spare Boss	Pack of 12	2,700
GLT-LW02-102	Laboratory Clamps — Size (Aprox.): For 1"(25mm) · Spare Boss	Pack of 6	1,950

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-09

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Filtration Assembly & Cannula S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

8
PRODUCTS

8
PACK
OPTIONS

₹4,800
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Filtration Assembly (S.S. 316 Quality)			
GLT-LW02-103	Filtration Assembly (S.S. 316 Quality) — FilterDia/Cup Capacity: 47mm/250ml	1 pc	10,750
GLT-LW02-104	Filtration Assembly (S.S. 316 Quality) — FilterDia/Cup Capacity: 47mm/500ml	1 pc	15,125
Cannula S.S.			
GLT-LW02-105	Cannula S.S. — Gauge: 16 · Dia - (ID): 1.6mm · Length: 18"	Pack of 12	4,800
GLT-LW02-106	Cannula S.S. — Gauge: 16 · Dia - (ID): 1.6mm · Length: 24"	Pack of 12	6,900
GLT-LW02-107	Cannula S.S. — Gauge: 16 · Dia - (ID): 1.6mm · Length: 30"	Pack of 12	9,000
GLT-LW02-108	Cannula S.S. — Gauge: 14 · Dia - (ID): 2mm · Length: 24"	Pack of 12	6,900
GLT-LW02-109	Cannula S.S. — Gauge: 14 · Dia - (ID): 2mm · Length: 30"	Pack of 12	9,000
GLT-LW02-110	Cannula S.S. — Gauge: 14 · Dia - (ID): 2mm · Length: 36"	Pack of 12	10,800

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-10

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Cannula S.S., Rack For Cryo Box S.S. & more

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

8
PRODUCTS

8
PACK
OPTIONS

₹5,250
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Cannula S.S.			
GLT-LW02-111	Cannula S.S. — Gauge: 12 · Dia - (ID): 2.5mm · Length: 24"	Pack of 12	6,900
GLT-LW02-112	Cannula S.S. — Gauge: 12 · Dia - (ID): 2.5mm · Length: 30"	Pack of 12	9,000
GLT-LW02-113	Cannula S.S. — Gauge: 12 · Dia - (ID): 2.5mm · Length: 36"	Pack of 12	10,800
Rack For Cryo Box S.S.			
GLT-LW02-114	Rack For Cryo Box S.S. — Number of Slot: 4 · Outer Size LxWx H(in mm): 145 × 145 × 250	1 pc	5,250
GLT-LW02-115	Rack For Cryo Box S.S. — Number of Slot: 6 · Outer Size LxWx H(in mm): 145 × 145 × 370	1 pc	5,750
GLT-LW02-116	Rack For Cryo Box S.S. — Number of Slot: 8 · Outer Size LxWx H(in mm): 145 × 145 × 500	1 pc	6,000
GLT-LW02-117	Rack For Cryo Box S.S. — Number of Slot: 10 · Outer Size LxWx H(in mm): 145 × 145 × 615	1 pc	6,500
Basket S.S.			
GLT-LW02-118	Basket S.S. — Number of Place: 4 Bottle · Capacity of Bottle: 500ml	1 pc	7,750

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-11

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Basket S.S., Tripod Stand S.S. & more

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

8
PRODUCTS

8
PACK
OPTIONS

₹1,200
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Basket S.S.			
GLT-LW02-119	Basket S.S. — Number of Place: 6 Bottle · Capacity of Bottle: 500ml	1 pc	8,250
GLT-LW02-120	Basket S.S. — Number of Place: 2 Bottle · Capacity of Bottle: 2500ml (2.5 Ltr)	1 pc	7,750
Tripod Stand S.S.			
GLT-LW02-121	Tripod Stand S.S. — Dia. x Height: 4" x 15" · Capacity of Flask For Separating Funnel: 1 Ltr	1 pc (sold in pack of 2)	1,500
GLT-LW02-122	Tripod Stand S.S. — Dia. x Height: 5" x 16" · Capacity of Flask For Separating Funnel: 2 Ltr	1 pc (sold in pack of 2)	1,875
GLT-LW02-123	Tripod Stand S.S. — Dia. x Height: 6" x 22" · Capacity of Flask For Separating Funnel: 3 Ltr	1 pc (sold in pack of 2)	2,375
GLT-LW02-124	Tripod Stand S.S. — Dia. x Height: 7" x 24" · Capacity of Flask For Separating Funnel: 4 Ltr	1 pc (sold in pack of 2)	2,750
Tripod Stand			
GLT-LW02-125	Tripod Stand — 6×4"	Pack of 6	1,200
GLT-LW02-126	Tripod Stand — 9×5"	Pack of 6	1,650

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-12

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Retort Ring & Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

10
PRODUCTS

10
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Retort Ring			
GLT-LW02-127	Retort Ring — Diameter of Glass: 2.5"	Pack of 6	1,200
GLT-LW02-128	Retort Ring — Diameter of Glass: 3"	Pack of 6	1,650
GLT-LW02-129	Retort Ring — Diameter of Glass: 4"	Pack of 6	2,325
GLT-LW02-130	Retort Ring — Diameter of Glass: 5"	Pack of 6	3,150
GLT-LW02-131	Retort Ring — Diameter of Glass: 6"	Pack of 6	3,600
Test Sieve			
GLT-LW02-132	Test Sieve — BSS 410/1969: 4 · MESH NO. ASTM 111/701: 5 · ISS 460/1972: 4.00mm · Sieve 4" dia.	1 pc	1,125
GLT-LW02-133	Test Sieve — BSS 410/1969: 4 · MESH NO. ASTM 111/701: 5 · ISS 460/1972: 4.00mm · Sieve 8" dia.	1 pc	1,375
GLT-LW02-134	Test Sieve — BSS 410/1969: 4 · MESH NO. ASTM 111/701: 5 · ISS 460/1972: 4.00mm · Sieve 12" dia.	1 pc	2,750
GLT-LW02-135	Test Sieve — BSS 410/1969: 4 · MESH NO. ASTM 111/701: 5 · ISS 460/1972: 4.00mm · Sieve 18" dia.	1 pc	4,500
GLT-LW02-136	Test Sieve — BSS 410/1969: 5 · MESH NO. ASTM 111/701: 6 · ISS 460/1972: 3.35mm · Sieve 4" dia.	1 pc	1,125

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-13

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

9
PRODUCTS

9
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-137	Test Sieve — BSS 410/1969: 5 · MESH NO. ASTM 111/701: 6 · ISS 460/1972: 3.35mm · Sieve 8" dia.	1 pc	1,375
GLT-LW02-138	Test Sieve — BSS 410/1969: 5 · MESH NO. ASTM 111/701: 6 · ISS 460/1972: 3.35mm · Sieve 12" dia.	1 pc	2,750
GLT-LW02-139	Test Sieve — BSS 410/1969: 5 · MESH NO. ASTM 111/701: 6 · ISS 460/1972: 3.35mm · Sieve 18" dia.	1 pc	4,500
GLT-LW02-140	Test Sieve — BSS 410/1969: 6 · MESH NO. ASTM 111/701: 7 · ISS 460/1972: 2.80mm · Sieve 4" dia.	1 pc	1,125
GLT-LW02-141	Test Sieve — BSS 410/1969: 6 · MESH NO. ASTM 111/701: 7 · ISS 460/1972: 2.80mm · Sieve 8" dia.	1 pc	1,375
GLT-LW02-142	Test Sieve — BSS 410/1969: 6 · MESH NO. ASTM 111/701: 7 · ISS 460/1972: 2.80mm · Sieve 12" dia.	1 pc	2,750
GLT-LW02-143	Test Sieve — BSS 410/1969: 6 · MESH NO. ASTM 111/701: 7 · ISS 460/1972: 2.80mm · Sieve 18" dia.	1 pc	4,500
GLT-LW02-144	Test Sieve — BSS 410/1969: 7 · MESH NO. ASTM 111/701: 8 · ISS 460/1972: 2.36mm · Sieve 4" dia.	1 pc	1,125
GLT-LW02-145	Test Sieve — BSS 410/1969: 7 · MESH NO. ASTM 111/701: 8 · ISS 460/1972: 2.36mm · Sieve 8" dia.	1 pc	1,375

Full specifications and every pack size: Product Catalogue 2026-27, pages 518-533.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-14

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

9
PRODUCTS

9
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-146	Test Sieve — BSS 410/1969: 7 · MESH NO. ASTM 111/701: 8 · ISS 460/1972: 2.36mm · Sieve 12" dia.	1 pc	2,750
GLT-LW02-147	Test Sieve — BSS 410/1969: 7 · MESH NO. ASTM 111/701: 8 · ISS 460/1972: 2.36mm · Sieve 18" dia.	1 pc	4,500
GLT-LW02-148	Test Sieve — BSS 410/1969: 8 · MESH NO. ASTM 111/701: 10 · ISS 460/1972: 2.00mm · Sieve 4" dia.	1 pc	1,125
GLT-LW02-149	Test Sieve — BSS 410/1969: 8 · MESH NO. ASTM 111/701: 10 · ISS 460/1972: 2.00mm · Sieve 8" dia.	1 pc	1,375
GLT-LW02-150	Test Sieve — BSS 410/1969: 8 · MESH NO. ASTM 111/701: 10 · ISS 460/1972: 2.00mm · Sieve 12" dia.	1 pc	2,750
GLT-LW02-151	Test Sieve — BSS 410/1969: 8 · MESH NO. ASTM 111/701: 10 · ISS 460/1972: 2.00mm · Sieve 18" dia.	1 pc	4,500
GLT-LW02-152	Test Sieve — BSS 410/1969: 10 · MESH NO. ASTM 111/701: 12 · ISS 460/1972: 1.70mm · Sieve 4" dia.	1 pc	1,125
GLT-LW02-153	Test Sieve — BSS 410/1969: 10 · MESH NO. ASTM 111/701: 12 · ISS 460/1972: 1.70mm · Sieve 8" dia.	1 pc	1,375
GLT-LW02-154	Test Sieve — BSS 410/1969: 10 · MESH NO. ASTM 111/701: 12 · ISS 460/1972: 1.70mm · Sieve 12" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-15

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

9
PRODUCTS

9
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-155	Test Sieve — BSS 410/1969: 10 · MESH NO. ASTM 111/701: 12 · ISS 460/1972: 1.70mm · Sieve 18" dia.	1 pc	4,500
GLT-LW02-156	Test Sieve — BSS 410/1969: 12 · MESH NO. ASTM 111/701: 14 · ISS 460/1972: 1.40mm · Sieve 4" dia.	1 pc	1,125
GLT-LW02-157	Test Sieve — BSS 410/1969: 12 · MESH NO. ASTM 111/701: 14 · ISS 460/1972: 1.40mm · Sieve 8" dia.	1 pc	1,375
GLT-LW02-158	Test Sieve — BSS 410/1969: 12 · MESH NO. ASTM 111/701: 14 · ISS 460/1972: 1.40mm · Sieve 12" dia.	1 pc	2,750
GLT-LW02-159	Test Sieve — BSS 410/1969: 12 · MESH NO. ASTM 111/701: 14 · ISS 460/1972: 1.40mm · Sieve 18" dia.	1 pc	4,500
GLT-LW02-160	Test Sieve — BSS 410/1969: 14 · MESH NO. ASTM 111/701: 16 · ISS 460/1972: 1.18mm · Sieve 4" dia.	1 pc	1,125
GLT-LW02-161	Test Sieve — BSS 410/1969: 14 · MESH NO. ASTM 111/701: 16 · ISS 460/1972: 1.18mm · Sieve 8" dia.	1 pc	1,375
GLT-LW02-162	Test Sieve — BSS 410/1969: 14 · MESH NO. ASTM 111/701: 16 · ISS 460/1972: 1.18mm · Sieve 12" dia.	1 pc	2,750
GLT-LW02-163	Test Sieve — BSS 410/1969: 14 · MESH NO. ASTM 111/701: 16 · ISS 460/1972: 1.18mm · Sieve 18" dia.	1 pc	4,500

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-16

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

8
PRODUCTS

8
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-164	Test Sieve — BSS 410/1969: 16 · MESH NO. ASTM 111/701: 18 · ISS 460/1972: 1.00mm · Sieve 4" dia.	1 pc	1,125
GLT-LW02-165	Test Sieve — BSS 410/1969: 16 · MESH NO. ASTM 111/701: 18 · ISS 460/1972: 1.00mm · Sieve 8" dia.	1 pc	1,375
GLT-LW02-166	Test Sieve — BSS 410/1969: 16 · MESH NO. ASTM 111/701: 18 · ISS 460/1972: 1.00mm · Sieve 12" dia.	1 pc	2,750
GLT-LW02-167	Test Sieve — BSS 410/1969: 16 · MESH NO. ASTM 111/701: 18 · ISS 460/1972: 1.00mm · Sieve 18" dia.	1 pc	4,500
GLT-LW02-168	Test Sieve — BSS 410/1969: 18 · MESH NO. ASTM 111/701: 20 · ISS 460/1972: 850 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-169	Test Sieve — BSS 410/1969: 18 · MESH NO. ASTM 111/701: 20 · ISS 460/1972: 850 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-170	Test Sieve — BSS 410/1969: 18 · MESH NO. ASTM 111/701: 20 · ISS 460/1972: 850 Micron · Sieve 12" dia.	1 pc	2,750
GLT-LW02-171	Test Sieve — BSS 410/1969: 18 · MESH NO. ASTM 111/701: 20 · ISS 460/1972: 850 Micron · Sieve 18" dia.	1 pc	4,500

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-17

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

7
PRODUCTS

7
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-172	Test Sieve — BSS 410/1969: 22 · MESH NO. ASTM 111/701: 25 · ISS 460/1972: 710 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-173	Test Sieve — BSS 410/1969: 22 · MESH NO. ASTM 111/701: 25 · ISS 460/1972: 710 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-174	Test Sieve — BSS 410/1969: 22 · MESH NO. ASTM 111/701: 25 · ISS 460/1972: 710 Micron · Sieve 12" dia.	1 pc	2,750
GLT-LW02-175	Test Sieve — BSS 410/1969: 22 · MESH NO. ASTM 111/701: 25 · ISS 460/1972: 710 Micron · Sieve 18" dia.	1 pc	4,500
GLT-LW02-176	Test Sieve — BSS 410/1969: 25 · MESH NO. ASTM 111/701: 30 · ISS 460/1972: 600 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-177	Test Sieve — BSS 410/1969: 25 · MESH NO. ASTM 111/701: 30 · ISS 460/1972: 600 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-178	Test Sieve — BSS 410/1969: 25 · MESH NO. ASTM 111/701: 30 · ISS 460/1972: 600 Micron · Sieve 12" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-18

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

7
PRODUCTS

7
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-179	Test Sieve — BSS 410/1969: 25 · MESH NO. ASTM 111/701: 30 · ISS 460/1972: 600 Micron · Sieve 18" dia.	1 pc	4,500
GLT-LW02-180	Test Sieve — BSS 410/1969: 30 · MESH NO. ASTM 111/701: 35 · ISS 460/1972: 500 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-181	Test Sieve — BSS 410/1969: 30 · MESH NO. ASTM 111/701: 35 · ISS 460/1972: 500 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-182	Test Sieve — BSS 410/1969: 30 · MESH NO. ASTM 111/701: 35 · ISS 460/1972: 500 Micron · Sieve 12" dia.	1 pc	2,750
GLT-LW02-183	Test Sieve — BSS 410/1969: 30 · MESH NO. ASTM 111/701: 35 · ISS 460/1972: 500 Micron · Sieve 18" dia.	1 pc	4,500
GLT-LW02-184	Test Sieve — BSS 410/1969: 36 · MESH NO. ASTM 111/701: 40 · ISS 460/1972: 425 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-185	Test Sieve — BSS 410/1969: 36 · MESH NO. ASTM 111/701: 40 · ISS 460/1972: 425 Micron · Sieve 8" dia.	1 pc	1,375

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-19

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

7
PRODUCTS

7
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-186	Test Sieve — BSS 410/1969: 36 · MESH NO. ASTM 111/701: 40 · ISS 460/1972: 425 Micron · Sieve 12" dia.	1 pc	2,750
GLT-LW02-187	Test Sieve — BSS 410/1969: 36 · MESH NO. ASTM 111/701: 40 · ISS 460/1972: 425 Micron · Sieve 18" dia.	1 pc	4,500
GLT-LW02-188	Test Sieve — BSS 410/1969: 44 · MESH NO. ASTM 111/701: 45 · ISS 460/1972: 355 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-189	Test Sieve — BSS 410/1969: 44 · MESH NO. ASTM 111/701: 45 · ISS 460/1972: 355 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-190	Test Sieve — BSS 410/1969: 44 · MESH NO. ASTM 111/701: 45 · ISS 460/1972: 355 Micron · Sieve 12" dia.	1 pc	2,750
GLT-LW02-191	Test Sieve — BSS 410/1969: 44 · MESH NO. ASTM 111/701: 45 · ISS 460/1972: 355 Micron · Sieve 18" dia.	1 pc	4,500
GLT-LW02-192	Test Sieve — BSS 410/1969: 52 · MESH NO. ASTM 111/701: 50 · ISS 460/1972: 300 Micron · Sieve 4" dia.	1 pc	1,125

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-20

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

7
PRODUCTS

7
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-193	Test Sieve — BSS 410/1969: 52 · MESH NO. ASTM 111/701: 50 · ISS 460/1972: 300 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-194	Test Sieve — BSS 410/1969: 52 · MESH NO. ASTM 111/701: 50 · ISS 460/1972: 300 Micron · Sieve 12" dia.	1 pc	2,750
GLT-LW02-195	Test Sieve — BSS 410/1969: 52 · MESH NO. ASTM 111/701: 50 · ISS 460/1972: 300 Micron · Sieve 18" dia.	1 pc	4,500
GLT-LW02-196	Test Sieve — BSS 410/1969: 60 · MESH NO. ASTM 111/701: 60 · ISS 460/1972: 250 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-197	Test Sieve — BSS 410/1969: 60 · MESH NO. ASTM 111/701: 60 · ISS 460/1972: 250 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-198	Test Sieve — BSS 410/1969: 60 · MESH NO. ASTM 111/701: 60 · ISS 460/1972: 250 Micron · Sieve 12" dia.	1 pc	2,750
GLT-LW02-199	Test Sieve — BSS 410/1969: 60 · MESH NO. ASTM 111/701: 60 · ISS 460/1972: 250 Micron · Sieve 18" dia.	1 pc	4,500

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-21

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

7
PRODUCTS

7
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-200	Test Sieve — BSS 410/1969: 72 · MESH NO. ASTM 111/701: 70 · ISS 460/1972: 212 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-201	Test Sieve — BSS 410/1969: 72 · MESH NO. ASTM 111/701: 70 · ISS 460/1972: 212 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-202	Test Sieve — BSS 410/1969: 72 · MESH NO. ASTM 111/701: 70 · ISS 460/1972: 212 Micron · Sieve 12" dia.	1 pc	2,750
GLT-LW02-203	Test Sieve — BSS 410/1969: 72 · MESH NO. ASTM 111/701: 70 · ISS 460/1972: 212 Micron · Sieve 18" dia.	1 pc	4,500
GLT-LW02-204	Test Sieve — BSS 410/1969: 85 · MESH NO. ASTM 111/701: 80 · ISS 460/1972: 180 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-205	Test Sieve — BSS 410/1969: 85 · MESH NO. ASTM 111/701: 80 · ISS 460/1972: 180 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-206	Test Sieve — BSS 410/1969: 85 · MESH NO. ASTM 111/701: 80 · ISS 460/1972: 180 Micron · Sieve 12" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-22

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹1,125
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-207	Test Sieve — BSS 410/1969: 85 · MESH NO. ASTM 111/701: 80 · ISS 460/1972: 180 Micron · Sieve 18" dia.	1 pc	4,500
GLT-LW02-208	Test Sieve — BSS 410/1969: 100 · MESH NO. ASTM 111/701: 100 · ISS 460/1972: 150 Micron · Sieve 4" dia.	1 pc	1,125
GLT-LW02-209	Test Sieve — BSS 410/1969: 100 · MESH NO. ASTM 111/701: 100 · ISS 460/1972: 150 Micron · Sieve 8" dia.	1 pc	1,375
GLT-LW02-210	Test Sieve — BSS 410/1969: 100 · MESH NO. ASTM 111/701: 100 · ISS 460/1972: 150 Micron · Sieve 12" dia.	1 pc	2,750
GLT-LW02-211	Test Sieve — BSS 410/1969: 100 · MESH NO. ASTM 111/701: 100 · ISS 460/1972: 150 Micron · Sieve 18" dia.	1 pc	4,500
GLT-LW02-212	Test Sieve — BSS 410/1969: 120 · MESH NO. ASTM 111/701: 120 · ISS 460/1972: 125 Micron · Sieve 4" dia.	1 pc	1,275

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-23

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹1,350
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-213	Test Sieve — BSS 410/1969: 120 · MESH NO. ASTM 111/701: 120 · ISS 460/1972: 125 Micron · Sieve 8" dia.	1 pc	1,600
GLT-LW02-214	Test Sieve — BSS 410/1969: 120 · MESH NO. ASTM 111/701: 120 · ISS 460/1972: 125 Micron · Sieve 12" dia.	1 pc	2,775
GLT-LW02-215	Test Sieve — BSS 410/1969: 120 · MESH NO. ASTM 111/701: 120 · ISS 460/1972: 125 Micron · Sieve 18" dia.	1 pc	4,825
GLT-LW02-216	Test Sieve — BSS 410/1969: 150 · MESH NO. ASTM 111/701: 140 · ISS 460/1972: 106 Micron · Sieve 4" dia.	1 pc	1,350
GLT-LW02-217	Test Sieve — BSS 410/1969: 150 · MESH NO. ASTM 111/701: 140 · ISS 460/1972: 106 Micron · Sieve 8" dia.	1 pc	1,700
GLT-LW02-218	Test Sieve — BSS 410/1969: 150 · MESH NO. ASTM 111/701: 140 · ISS 460/1972: 106 Micron · Sieve 12" dia.	1 pc	2,775

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-24

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹1,512
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-219	Test Sieve — BSS 410/1969: 150 · MESH NO. ASTM 111/701: 140 · ISS 460/1972: 106 Micron · Sieve 18" dia.	1 pc	4,825
GLT-LW02-220	Test Sieve — BSS 410/1969: 170 · MESH NO. ASTM 111/701: 170 · ISS 460/1972: 90 Micron · Sieve 4" dia.	1 pc	1,512
GLT-LW02-221	Test Sieve — BSS 410/1969: 170 · MESH NO. ASTM 111/701: 170 · ISS 460/1972: 90 Micron · Sieve 8" dia.	1 pc	1,900
GLT-LW02-222	Test Sieve — BSS 410/1969: 170 · MESH NO. ASTM 111/701: 170 · ISS 460/1972: 90 Micron · Sieve 12" dia.	1 pc	3,075
GLT-LW02-223	Test Sieve — BSS 410/1969: 170 · MESH NO. ASTM 111/701: 170 · ISS 460/1972: 90 Micron · Sieve 18" dia.	1 pc	5,075
GLT-LW02-224	Test Sieve — BSS 410/1969: 200 · MESH NO. ASTM 111/701: 200 · ISS 460/1972: 75 Micron · Sieve 4" dia.	1 pc	1,950

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-25

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,025
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-225	Test Sieve — BSS 410/1969: 200 · MESH NO. ASTM 111/701: 200 · ISS 460/1972: 75 Micron · Sieve 8" dia.	1 pc	2,325
GLT-LW02-226	Test Sieve — BSS 410/1969: 200 · MESH NO. ASTM 111/701: 200 · ISS 460/1972: 75 Micron · Sieve 12" dia.	1 pc	3,150
GLT-LW02-227	Test Sieve — BSS 410/1969: 200 · MESH NO. ASTM 111/701: 200 · ISS 460/1972: 75 Micron · Sieve 18" dia.	1 pc	5,275
GLT-LW02-228	Test Sieve — BSS 410/1969: 240 · MESH NO. ASTM 111/701: 230 · ISS 460/1972: 63 Micron · Sieve 4" dia.	1 pc	2,025
GLT-LW02-229	Test Sieve — BSS 410/1969: 240 · MESH NO. ASTM 111/701: 230 · ISS 460/1972: 63 Micron · Sieve 8" dia.	1 pc	2,650
GLT-LW02-230	Test Sieve — BSS 410/1969: 240 · MESH NO. ASTM 111/701: 230 · ISS 460/1972: 63 Micron · Sieve 12" dia.	1 pc	3,300

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-26

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,325
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-231	Test Sieve — BSS 410/1969: 240 · MESH NO. ASTM 111/701: 230 · ISS 460/1972: 63 Micron · Sieve 18" dia.	1 pc	5,575
GLT-LW02-232	Test Sieve — BSS 410/1969: 300 · MESH NO. ASTM 111/701: 270 · ISS 460/1972: 53 Micron · Sieve 4" dia.	1 pc	2,325
GLT-LW02-233	Test Sieve — BSS 410/1969: 300 · MESH NO. ASTM 111/701: 270 · ISS 460/1972: 53 Micron · Sieve 8" dia.	1 pc	3,075
GLT-LW02-234	Test Sieve — BSS 410/1969: 300 · MESH NO. ASTM 111/701: 270 · ISS 460/1972: 53 Micron · Sieve 12" dia.	1 pc	3,525
GLT-LW02-235	Test Sieve — BSS 410/1969: 300 · MESH NO. ASTM 111/701: 270 · ISS 460/1972: 53 Micron · Sieve 18" dia.	1 pc	5,625
GLT-LW02-236	Test Sieve — BSS 410/1969: 350 · MESH NO. ASTM 111/701: 325 · ISS 460/1972: 45 Micron · 18"DIA: — Sieve 4" dia.	1 pc	2,725

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-27

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹3,012
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-237	Test Sieve — BSS 410/1969: 350 · MESH NO. ASTM 111/701: 325 · ISS 460/1972: 45 Micron · 18"DIA: — · Sieve 8" dia.	1 pc	3,512
GLT-LW02-238	Test Sieve — BSS 410/1969: 350 · MESH NO. ASTM 111/701: 325 · ISS 460/1972: 45 Micron · 18"DIA: — · Sieve 12" dia.	1 pc	6,500
GLT-LW02-239	Test Sieve — BSS 410/1969: 400 · MESH NO. ASTM 111/701: 400 · ISS 460/1972: 37 Micron · 18"DIA: — · Sieve 4" dia.	1 pc	3,012
GLT-LW02-240	Test Sieve — BSS 410/1969: 400 · MESH NO. ASTM 111/701: 400 · ISS 460/1972: 37 Micron · 18"DIA: — · Sieve 8" dia.	1 pc	4,075
GLT-LW02-241	Test Sieve — BSS 410/1969: 400 · MESH NO. ASTM 111/701: 400 · ISS 460/1972: 37 Micron · 18"DIA: — · Sieve 12" dia.	1 pc	7,000
GLT-LW02-242	Test Sieve — BSS 410/1969: 500 · MESH NO. ASTM 111/701: — · ISS 460/1972: 25 Micron · 18"DIA: — · Sieve 4" dia.	1 pc	4,600

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-28

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹1,500
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve			
GLT-LW02-243	Test Sieve — BSS 410/1969: 500 · MESH NO. ASTM 111/701: - · ISS 460/1972: 25 Micron · 18"DIA: - · Sieve 8" dia.	1 pc	6,500
GLT-LW02-244	Test Sieve — BSS 410/1969: 500 · MESH NO. ASTM 111/701: - · ISS 460/1972: 25 Micron · 18"DIA: - · Sieve 12" dia.	1 pc	On request
GLT-LW02-245	Test Sieve — BSS 410/1969: Lid · MESH NO. ASTM 111/701: &Receiver For · ISS 460/1972: Above Per Set · Sieve 4" dia.	1 pc	1,500
GLT-LW02-246	Test Sieve — BSS 410/1969: Lid · MESH NO. ASTM 111/701: &Receiver For · ISS 460/1972: Above Per Set · Sieve 8" dia.	1 pc	1,562
GLT-LW02-247	Test Sieve — BSS 410/1969: Lid · MESH NO. ASTM 111/701: &Receiver For · ISS 460/1972: Above Per Set · Sieve 12" dia.	1 pc	1,950
GLT-LW02-248	Test Sieve — BSS 410/1969: Lid · MESH NO. ASTM 111/701: &Receiver For · ISS 460/1972: Above Per Set · Sieve 18" dia.	1 pc	4,125

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-29

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

9
PRODUCTS

9
PACK
OPTIONS

₹2,750
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-249	Test Sieve S.S. — BSS 410/1969: 4 · MESH NO. ASTM 111/701: 5 · ISS 460/1972: 4.00mm · Sieve 8" dia.	1 pc	2,750
GLT-LW02-250	Test Sieve S.S. — BSS 410/1969: 4 · MESH NO. ASTM 111/701: 5 · ISS 460/1972: 4.00mm · Sieve 12" dia.	1 pc	4,500
GLT-LW02-251	Test Sieve S.S. — BSS 410/1969: 4 · MESH NO. ASTM 111/701: 5 · ISS 460/1972: 4.00mm · Sieve 18" dia.	1 pc	6,000
GLT-LW02-252	Test Sieve S.S. — BSS 410/1969: 5 · MESH NO. ASTM 111/701: 6 · ISS 460/1972: 3.35mm · Sieve 8" dia.	1 pc	2,750
GLT-LW02-253	Test Sieve S.S. — BSS 410/1969: 5 · MESH NO. ASTM 111/701: 6 · ISS 460/1972: 3.35mm · Sieve 12" dia.	1 pc	4,500
GLT-LW02-254	Test Sieve S.S. — BSS 410/1969: 5 · MESH NO. ASTM 111/701: 6 · ISS 460/1972: 3.35mm · Sieve 18" dia.	1 pc	6,000
GLT-LW02-255	Test Sieve S.S. — BSS 410/1969: 6 · MESH NO. ASTM 111/701: 7 · ISS 460/1972: 2.80mm · Sieve 8" dia.	1 pc	2,750
GLT-LW02-256	Test Sieve S.S. — BSS 410/1969: 6 · MESH NO. ASTM 111/701: 7 · ISS 460/1972: 2.80mm · Sieve 12" dia.	1 pc	4,500
GLT-LW02-257	Test Sieve S.S. — BSS 410/1969: 6 · MESH NO. ASTM 111/701: 7 · ISS 460/1972: 2.80mm · Sieve 18" dia.	1 pc	6,000

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-30

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

7
PRODUCTS

7
PACK
OPTIONS

₹2,750
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-258	Test Sieve S.S. — BSS 410/1969: 7 · MESH NO. ASTM 111/701: 8 · ISS 460/1972: 2.36mm · Sieve 8" dia.	1 pc	2,750
GLT-LW02-259	Test Sieve S.S. — BSS 410/1969: 7 · MESH NO. ASTM 111/701: 8 · ISS 460/1972: 2.36mm · Sieve 12" dia.	1 pc	4,500
GLT-LW02-260	Test Sieve S.S. — BSS 410/1969: 7 · MESH NO. ASTM 111/701: 8 · ISS 460/1972: 2.36mm · Sieve 18" dia.	1 pc	6,000
GLT-LW02-261	Test Sieve S.S. — BSS 410/1969: 8 · MESH NO. ASTM 111/701: 10 · ISS 460/1972: 2.00mm · Sieve 8" dia.	1 pc	2,750
GLT-LW02-262	Test Sieve S.S. — BSS 410/1969: 8 · MESH NO. ASTM 111/701: 10 · ISS 460/1972: 2.00mm · Sieve 12" dia.	1 pc	4,500
GLT-LW02-263	Test Sieve S.S. — BSS 410/1969: 8 · MESH NO. ASTM 111/701: 10 · ISS 460/1972: 2.00mm · Sieve 18" dia.	1 pc	6,000
GLT-LW02-264	Test Sieve S.S. — BSS 410/1969: 10 · MESH NO. ASTM 111/701: 12 · ISS 460/1972: 1.70mm · Sieve 8" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-31

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,750
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-265	Test Sieve S.S. — BSS 410/1969: 10 · MESH NO. ASTM 111/701: 12 · ISS 460/1972: 1.70mm · Sieve 12" dia.	1 pc	4,500
GLT-LW02-266	Test Sieve S.S. — BSS 410/1969: 10 · MESH NO. ASTM 111/701: 12 · ISS 460/1972: 1.70mm · Sieve 18" dia.	1 pc	6,000
GLT-LW02-267	Test Sieve S.S. — BSS 410/1969: 12 · MESH NO. ASTM 111/701: 14 · ISS 460/1972: 1.40mm · Sieve 8" dia.	1 pc	2,750
GLT-LW02-268	Test Sieve S.S. — BSS 410/1969: 12 · MESH NO. ASTM 111/701: 14 · ISS 460/1972: 1.40mm · Sieve 12" dia.	1 pc	4,500
GLT-LW02-269	Test Sieve S.S. — BSS 410/1969: 12 · MESH NO. ASTM 111/701: 14 · ISS 460/1972: 1.40mm · Sieve 18" dia.	1 pc	6,000
GLT-LW02-270	Test Sieve S.S. — BSS 410/1969: 14 · MESH NO. ASTM 111/701: 16 · ISS 460/1972: 1.18mm · Sieve 8" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-32

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,750
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-271	Test Sieve S.S. — BSS 410/1969: 14 · MESH NO. ASTM 111/701: 16 · ISS 460/1972: 1.18mm · Sieve 12" dia.	1 pc	4,500
GLT-LW02-272	Test Sieve S.S. — BSS 410/1969: 14 · MESH NO. ASTM 111/701: 16 · ISS 460/1972: 1.18mm · Sieve 18" dia.	1 pc	6,000
GLT-LW02-273	Test Sieve S.S. — BSS 410/1969: 16 · MESH NO. ASTM 111/701: 18 · ISS 460/1972: 1.00mm · Sieve 8" dia.	1 pc	2,750
GLT-LW02-274	Test Sieve S.S. — BSS 410/1969: 16 · MESH NO. ASTM 111/701: 18 · ISS 460/1972: 1.00mm · Sieve 12" dia.	1 pc	4,500
GLT-LW02-275	Test Sieve S.S. — BSS 410/1969: 16 · MESH NO. ASTM 111/701: 18 · ISS 460/1972: 1.00mm · Sieve 18" dia.	1 pc	6,000
GLT-LW02-276	Test Sieve S.S. — BSS 410/1969: 18 · MESH NO. ASTM 111/701: 20 · ISS 460/1972: 850 Micro · Sieve 8" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-33

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,750
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-277	Test Sieve S.S. — BSS 410/1969: 18 · MESH NO. ASTM 111/701: 20 · ISS 460/1972: 850 Micro · Sieve 12" dia.	1 pc	4,500
GLT-LW02-278	Test Sieve S.S. — BSS 410/1969: 18 · MESH NO. ASTM 111/701: 20 · ISS 460/1972: 850 Micro · Sieve 18" dia.	1 pc	6,000
GLT-LW02-279	Test Sieve S.S. — BSS 410/1969: 22 · MESH NO. ASTM 111/701: 25 · ISS 460/1972: 710 Micron · Sieve 8" dia.	1 pc	2,750
GLT-LW02-280	Test Sieve S.S. — BSS 410/1969: 22 · MESH NO. ASTM 111/701: 25 · ISS 460/1972: 710 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-281	Test Sieve S.S. — BSS 410/1969: 22 · MESH NO. ASTM 111/701: 25 · ISS 460/1972: 710 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-282	Test Sieve S.S. — BSS 410/1969: 25 · MESH NO. ASTM 111/701: 30 · ISS 460/1972: 600 Micron · Sieve 8" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-34

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,750
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-283	Test Sieve S.S. — BSS 410/1969: 25 · MESH NO. ASTM 111/701: 30 · ISS 460/1972: 600 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-284	Test Sieve S.S. — BSS 410/1969: 25 · MESH NO. ASTM 111/701: 30 · ISS 460/1972: 600 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-285	Test Sieve S.S. — BSS 410/1969: 30 · MESH NO. ASTM 111/701: 35 · ISS 460/1972: 500 Micron · Sieve 8" dia.	1 pc	2,750
GLT-LW02-286	Test Sieve S.S. — BSS 410/1969: 30 · MESH NO. ASTM 111/701: 35 · ISS 460/1972: 500 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-287	Test Sieve S.S. — BSS 410/1969: 30 · MESH NO. ASTM 111/701: 35 · ISS 460/1972: 500 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-288	Test Sieve S.S. — BSS 410/1969: 36 · MESH NO. ASTM 111/701: 40 · ISS 460/1972: 425 Micron · Sieve 8" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-35

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,750
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-289	Test Sieve S.S. — BSS 410/1969: 36 · MESH NO. ASTM 111/701: 40 · ISS 460/1972: 425 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-290	Test Sieve S.S. — BSS 410/1969: 36 · MESH NO. ASTM 111/701: 40 · ISS 460/1972: 425 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-291	Test Sieve S.S. — BSS 410/1969: 44 · MESH NO. ASTM 111/701: 45 · ISS 460/1972: 355 Micron · Sieve 8" dia.	1 pc	2,750
GLT-LW02-292	Test Sieve S.S. — BSS 410/1969: 44 · MESH NO. ASTM 111/701: 45 · ISS 460/1972: 355 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-293	Test Sieve S.S. — BSS 410/1969: 44 · MESH NO. ASTM 111/701: 45 · ISS 460/1972: 355 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-294	Test Sieve S.S. — BSS 410/1969: 52 · MESH NO. ASTM 111/701: 50 · ISS 460/1972: 300 Micron · Sieve 8" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-36

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,750
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-295	Test Sieve S.S. — BSS 410/1969: 52 · MESH NO. ASTM 111/701: 50 · ISS 460/1972: 300 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-296	Test Sieve S.S. — BSS 410/1969: 52 · MESH NO. ASTM 111/701: 50 · ISS 460/1972: 300 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-297	Test Sieve S.S. — BSS 410/1969: 60 · MESH NO. ASTM 111/701: 60 · ISS 460/1972: 250 Micron · Sieve 8" dia.	1 pc	2,750
GLT-LW02-298	Test Sieve S.S. — BSS 410/1969: 60 · MESH NO. ASTM 111/701: 60 · ISS 460/1972: 250 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-299	Test Sieve S.S. — BSS 410/1969: 60 · MESH NO. ASTM 111/701: 60 · ISS 460/1972: 250 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-300	Test Sieve S.S. — BSS 410/1969: 72 · MESH NO. ASTM 111/701: 70 · ISS 460/1972: 212 Micron · Sieve 8" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-37

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,750
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-301	Test Sieve S.S. — BSS 410/1969: 72 · MESH NO. ASTM 111/701: 70 · ISS 460/1972: 212 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-302	Test Sieve S.S. — BSS 410/1969: 72 · MESH NO. ASTM 111/701: 70 · ISS 460/1972: 212 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-303	Test Sieve S.S. — BSS 410/1969: 85 · MESH NO. ASTM 111/701: 80 · ISS 460/1972: 180 Micron · Sieve 8" dia.	1 pc	2,750
GLT-LW02-304	Test Sieve S.S. — BSS 410/1969: 85 · MESH NO. ASTM 111/701: 80 · ISS 460/1972: 180 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-305	Test Sieve S.S. — BSS 410/1969: 85 · MESH NO. ASTM 111/701: 80 · ISS 460/1972: 180 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-306	Test Sieve S.S. — BSS 410/1969: 100 · MESH NO. ASTM 111/701: 100 · ISS 460/1972: 150 Micron · Sieve 8" dia.	1 pc	2,750

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-38

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹2,938
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-307	Test Sieve S.S. — BSS 410/1969: 100 · MESH NO. ASTM 111/701: 100 · ISS 460/1972: 150 Micron · Sieve 12" dia.	1 pc	4,500
GLT-LW02-308	Test Sieve S.S. — BSS 410/1969: 100 · MESH NO. ASTM 111/701: 100 · ISS 460/1972: 150 Micron · Sieve 18" dia.	1 pc	6,000
GLT-LW02-309	Test Sieve S.S. — BSS 410/1969: 120 · MESH NO. ASTM 111/701: 120 · ISS 460/1972: 125 Micron · Sieve 8" dia.	1 pc	2,938
GLT-LW02-310	Test Sieve S.S. — BSS 410/1969: 120 · MESH NO. ASTM 111/701: 120 · ISS 460/1972: 125 Micron · Sieve 12" dia.	1 pc	4,688
GLT-LW02-311	Test Sieve S.S. — BSS 410/1969: 120 · MESH NO. ASTM 111/701: 120 · ISS 460/1972: 125 Micron · Sieve 18" dia.	1 pc	6,188
GLT-LW02-312	Test Sieve S.S. — BSS 410/1969: 150 · MESH NO. ASTM 111/701: 140 · ISS 460/1972: 106 Micron · Sieve 8" dia.	1 pc	3,100

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-39

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹3,325
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-313	Test Sieve S.S. — BSS 410/1969: 150 · MESH NO. ASTM 111/701: 140 · ISS 460/1972: 106 Micron · Sieve 12" dia.	1 pc	4,825
GLT-LW02-314	Test Sieve S.S. — BSS 410/1969: 150 · MESH NO. ASTM 111/701: 140 · ISS 460/1972: 106 Micron · Sieve 18" dia.	1 pc	6,575
GLT-LW02-315	Test Sieve S.S. — BSS 410/1969: 170 · MESH NO. ASTM 111/701: 170 · ISS 460/1972: 90 Micron · Sieve 8" dia.	1 pc	3,325
GLT-LW02-316	Test Sieve S.S. — BSS 410/1969: 170 · MESH NO. ASTM 111/701: 170 · ISS 460/1972: 90 Micron · Sieve 12" dia.	1 pc	5,025
GLT-LW02-317	Test Sieve S.S. — BSS 410/1969: 170 · MESH NO. ASTM 111/701: 170 · ISS 460/1972: 90 Micron · Sieve 18" dia.	1 pc	7,025
GLT-LW02-318	Test Sieve S.S. — BSS 410/1969: 200 · MESH NO. ASTM 111/701: 200 · ISS 460/1972: 75 Micron · Sieve 8" dia.	1 pc	3,525

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-40

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹3,825
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-319	Test Sieve S.S. — BSS 410/1969: 200 · MESH NO. ASTM 111/701: 200 · ISS 460/1972: 75 Micron · Sieve 12" dia.	1 pc	5,250
GLT-LW02-320	Test Sieve S.S. — BSS 410/1969: 200 · MESH NO. ASTM 111/701: 200 · ISS 460/1972: 75 Micron · Sieve 18" dia.	1 pc	7,275
GLT-LW02-321	Test Sieve S.S. — BSS 410/1969: 240 · MESH NO. ASTM 111/701: 230 · ISS 460/1972: 63 Micron · Sieve 8" dia.	1 pc	3,825
GLT-LW02-322	Test Sieve S.S. — BSS 410/1969: 240 · MESH NO. ASTM 111/701: 230 · ISS 460/1972: 63 Micron · Sieve 12" dia.	1 pc	5,575
GLT-LW02-323	Test Sieve S.S. — BSS 410/1969: 240 · MESH NO. ASTM 111/701: 230 · ISS 460/1972: 63 Micron · Sieve 18" dia.	1 pc	7,750
GLT-LW02-324	Test Sieve S.S. — BSS 410/1969: 300 · MESH NO. ASTM 111/701: 270 · ISS 460/1972: 53 Micron · Sieve 8" dia.	1 pc	4,525

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-41

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S.

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

6
PRODUCTS

6
PACK
OPTIONS

₹4,950
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-325	Test Sieve S.S. — BSS 410/1969: 300 · MESH NO. ASTM 111/701: 270 · ISS 460/1972: 53 Micron · Sieve 12" dia.	1 pc	5,850
GLT-LW02-326	Test Sieve S.S. — BSS 410/1969: 300 · MESH NO. ASTM 111/701: 270 · ISS 460/1972: 53 Micron · Sieve 18" dia.	1 pc	8,500
GLT-LW02-327	Test Sieve S.S. — BSS 410/1969: 350 · MESH NO. ASTM 111/701: 325 · ISS 460/1972: 45 Micron · 18"DIA: - · Sieve 8" dia.	1 pc	4,950
GLT-LW02-328	Test Sieve S.S. — BSS 410/1969: 350 · MESH NO. ASTM 111/701: 325 · ISS 460/1972: 45 Micron · 18"DIA: - · Sieve 12" dia.	1 pc	7,375
GLT-LW02-329	Test Sieve S.S. — BSS 410/1969: 400 · MESH NO. ASTM 111/701: 400 · ISS 460/1972: 37 Micron · 18"DIA: - · Sieve 8" dia.	1 pc	5,500
GLT-LW02-330	Test Sieve S.S. — BSS 410/1969: 400 · MESH NO. ASTM 111/701: 400 · ISS 460/1972: 37 Micron · 18"DIA: - · Sieve 12" dia.	1 pc	8,000

Full specifications and every pack size: Product Catalogue 2026-27, pages **518-533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026-27 · Edition 2.0 ·

GLT-FLY-LW02-42

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113

LW02 STAINLESS STEEL & METAL LABWARE

Test Sieve S.S., Bench Top Beta Shield & more

Illustrative image

Stainless-steel trays, stands, racks, beakers, instruments, clamps, lab jacks and test sieves.

8
PRODUCTS

8
PACK
OPTIONS

₹325
PRICE FROM

TYPICAL RESEARCH USE

Holding and supporting apparatus

Dissection and sample handling

Particle sizing with sieves

GRACE CAT. NO.	PRODUCT	PACK SIZES	PRICE (₹)
Test Sieve S.S.			
GLT-LW02-331	Test Sieve S.S. — BSS 410/1969: 500 · MESH NO. ASTM 111/701: — · ISS 460/1972: 25 Micron · 18"DIA: — · Sieve 8" dia.	1 pc	9,250
GLT-LW02-332	Test Sieve S.S. — BSS 410/1969: 500 · MESH NO. ASTM 111/701: — · ISS 460/1972: 25 Micron · 18"DIA: — · Sieve 12" dia.	1 pc	On request
Bench Top Beta Shield			
GLT-LW02-333	Bench Top Beta Shield — 45 × 30 × 30	1 pc	10,738
Plastic Forceps			
GLT-LW02-334	Plastic Forceps — Plastic forceps · 5"	Pack of 12	325
Spatula Plastic			
GLT-LW02-335	Spatula Plastic — 150mm	Pack of 12	450
GLT-LW02-336	Spatula Plastic — 200mm	Pack of 12	675
Staining Box			
GLT-LW02-337	Staining Box — Size in mm(1 xwxh): 125×125×50	Pack of 2	1,400
GLT-LW02-338	Staining Box — Size in mm(1 xwxh): 225×225×50	Pack of 2	2,550

Full specifications and every pack size: Product Catalogue 2026–27, pages **518–533**.

Grace Lifetech Pvt. Ltd.

Reg. Office: Mehrauli, New Delhi, India · www.gracelifetech.com

Order by Grace Cat. No. · quotations on request

Prices in INR, exclusive of GST and freight; subject to change without notice. "From" = lowest listed pack price.

Product Catalogue 2026–27 · Edition 2.0 ·

GLT-FLY-LW02-43

R&D Centre: TIH/CUB/04, BTCIF, 1st Floor, TICEL Biopark Limited Phase II, No. 5, CSIR Road, Taramani, Chennai – 600113